


COMPO-VIBES

Smile Line 
powered by : **STYLEITALIANO**

Smart vibrations !

The guarantee of an easier, faster and better composite's modeling !
Thanks to the micro-vibrations it generates, Compo-Vibes optimises the material's thixotropic properties and makes your work easier :


- Better adaptation of the composite to the cavity wall, no detachment, less air bubbles
- Composite much less sticky than with a traditional instrument
- Modeling and sculpture of composites much easier and fast
- Increased control of material thicknesses
- More precise application



Watch the video !



COMPO-VIBES



Technical data :

- Compact and light 160 x 13mm / Weight 33g
- Button allowing a vibration in continuous mode or with short impulses
- Four instruments, very easily exchangeable
- No wire, AAAA battery 1.5V
- Vibration frequency 140 Hz

Description	Item no.
Compo-Vibes Basic Set : Compo-Vibes unit with one instrument «Tip medium»	03000-CV-Basic
Compo-Vibes Master Set : Compo-Vibes unit with 4 instruments «Spatula / Tip large / Tip medium / tip thin »	03000-CV-Master
Set of 4 instruments : « Spatula / Tip large / Tip medium / Tip thin »	03010-CV
Spatula, 2 pcs	03011-CV
Instrument "Tip large", 2 pcs	03012-CV
Instrument "Tip medium", 2 pcs	03013-CV
Instrument "Tip thin", 2 pcs	03014-CV
Battery AAAA 1,5 V	01030-MV
Disposable protective sleeves, 200 pcs	3011-KL

i create... we create