

Benefits of upgrading to the all NEW STRATA-G™



All-New Separator Rings

- **Lasts longer!** Drawn-wire nickel-titanium holds its shape better.
- **Enhanced marginal ridge!** Helps support the matrix band in the proper position.
- **Smaller size!** Stays out of your way.




Wide Preparation Ring

- **Garrison exclusive!** Only Garrison helps you restore the largest preps.


Full Curve Matrix Bands

- **Enhanced marginal ridge!** Helps form this critical anatomy.
- **Improved emergence profile!** Ideal anatomy is built right in.
- **Subgingival extensions!** All four of the standard sizes are ready for deeper preps.
- **Still dead-soft** and only 0.0016" thick.



Ultra-Adaptive Wedges

- **Smooth tunnel design!** Sits lower interproximally and seals deeper.
- **More flexibility!** The Soft-Face™ over-mold with fins and flexible sidewalls easily adapt to concavities and irregularities.
- **Stackable!** Secure stacking of multiple wedges when simultaneous wedging from the buccal and lingual is desired.



Dentistry's Premier Forceps

- **Larger offset!** Makes ring placement easier.
- **Double tipped!** Place rings in either direction for maximum access.
- **Places virtually any ring from any manufacturer!**



SGR-KSH-10 contains:

- 1 Short Ring (blue)
- 1 Tall Ring (orange)
- 1 Wide Prep Ring (green)
- 50 Assorted Matrix Bands
- 50 Assorted Firm Matrix Bands
- 10 VariStrip™
- 80 Assorted Wedges
- 12 Assorted Rally™ Polishers
- 1 Ring Placement Forceps (FXP01)

ORDER TODAY! 888.437.0032 or contact your authorized Garrison dealer.

Garrison

Distributed by Curion
curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200

CURION
dentistry. curated.