



Composite Luting Cement

DUO-LINK is a dual-cured luting cement with ultrafine filler particles and is intended to be used with an adhesive system. Available in two shades: TRANSLUCENT and CLEAR, the exceptionally high glass filler content increases strength making DUO-LINK resistant to wear and staining. It is ideal for cementing all-ceramic/porcelain, metal, and composite restorations (inlays, onlays, crowns), as well as fiber and metal posts.

- **1 DUAL-SYRINGE PACKAGE** (Translucent or Clear)
 1 Dual-Syringe (Base/Catalyst, 8g mixed), 15 Mixing Tips and Intraoral Tips, Instructions/MSDS
- **2 SYRINGE PACKAGE** (Translucent)
 1 Syringe Base (4g), 1 Syringe Catalyst (4g), Accessories, Instructions/MSDS



UNIQUE BENEFITS

- Extremely high degree of conversion in both light-cured and self-cured modes ensuring maximum performance for any clinical situation
- Exceptional high glass filler content increases strength making DUO-LINK resistant to wear
- Low film thickness ensures the restoration is completely seated
- Auto-mix dual-syringe guarantees a perfect mix, easy placement and limited waste

CLINICAL SIGNIFICANCE

- Extremely high degree of conversion and stability ensures a long lasting restoration placement
- A universal cement with a proven clinical track record

PHYSICAL PROPERTIES

% Filled (Weight/Volume)	Particle Size	Film Thickness	Diametral Tensile Strength	Compressive Strength	Flexural Strength
66/42	<1.0µm	17µm	49 MPa	228 MPa	126 MPa

+ RELATED MATERIALS



Double Taper LIGHT-POST™ ILLUSION™



ONE-STEP® / ONE-STEP® PLUS



BISCO Etchants