

Composi-Tight® 3D Fusion™

Firm Matrix Bands

NEW!

New **Composi-Tight® 3D Fusion™** Firm matrix bands with **eZ-Place™** technology maintain their shape during placement making them *significantly easier* to insert!

eZ-Place™ technology combines specially processed stainless steel with a proprietary matrix band forming process to yield Fusion Firm's remarkable handling characteristics. It resists deformation nearly 70% better than typical dead-soft bands yet is pre-contoured and super thin at 0.0015". You're ensured the tight contacts and ideal anatomy you've come to expect from the Composi-Tight® 3D Fusion™ system.

Fusion Firm:

- Nearly 70% better at resisting deformation!
- Super-thin at only 0.0015"
- Five ideal sizes!

Perfect for:

- Tight interproximal spaces!
- More conservative preps!
- First-time users and veterans alike!

Composi-Tight® 3D Fusion™ Firm matrix bands are now available in a Composi-Tight® 3D Fusion™ kit.

FXR-HHF-10

Composi-Tight® 3D Fusion™ Firm Bands

Description

Assorted Band Kit: 5 sizes

Assorted Band Mini Kit: 5 sizes

Item

FXHB04

FXHB05

| FXH Series | Composi-Tight® 3D Fusion™ Firm | | Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement. | | | | | Item |
|------------|--------------------------------|--------|---|------------------|----------------------------------|------------------|------------------|----------------------------------|
| | thickness 0.0015" | height | NA | | | | | |
| | | | | FXH100 4.3 mm | FXH150 4.3 mm w/ext 5.9 mm | FXH175 5.3 mm | FXH200 6.1 mm | FXH300 6.1 mm w/ext 7.7 mm |

Item

FXH100

FXH150

FXH175

FXH200

FXH300

Be sure to check out the website at www.garrisondental.com for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.


Distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200




























curion
dentistry. curated.

Garrison BAND COMPARISON CHART

Dental Solutions



FX Series
FXH Series
SXR Series
B Series
M - SM - AU Series

| | | | | | | |
|------------|--------------------------------|---|--|--|--|--|
| FXH Series | Composi-Tight® 3D Fusion™ Firm | Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement. | | | | |
| | thickness 0.0015" | NA |  FXH100 4.3 mm |  FXH150 4.3 mm w/ext 5.9 mm |  FXH175 5.3 mm |  FXH200 6.1 mm |
| FX Series | Composi-Tight® 3D Fusion™ | Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab. | | | | |
| | thickness 0.0016" | NA |  FX100 4.0 mm w/ext 4.4 mm |  FX150 4.1 mm w/ext 6.0 mm |  FX175 5.0 mm w/ext 5.6 mm |  FX200 5.4 mm w/ext 6.6 mm |
| SXR Series | Slick Bands™ XR | Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement. | | | | |
| | thickness 0.0016" | NA |  SXR100 4.6 mm |  SXR150 3.8 mm w/ext 5.4 mm |  SXR175 5.4 mm |  SXR200 6.4 mm |
| SM Series | Slick Bands™ | Slick Bands™ Matrices have non-stick coating. | | | | |
| | thickness 0.0016" |  SM050 3.2 mm w/ext 4.2 mm |  SM100 4.6 mm |  SM150 3.8 mm w/ext 5.4 mm |  SM175 5.5 mm |  SM200 6.4 mm |
| M Series | Composi-Tight® | Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour. | | | | |
| | thickness 0.0015" |  M050 3.2 mm w/ext 4.2 mm |  M100 4.6 mm |  M150 3.8 mm w/ext 5.4 mm |  M175 5.5 mm |  M200 6.4 mm |
| AU Series | Composi-Tight® Gold® | The Gold bands are the same as the Composi-Tight® bands. | | | | |
| | thickness 0.0015" |  AU050 3.2 mm w/ext 4.2 mm |  AU100 4.6 mm |  AU150 3.8 mm w/ext 5.4 mm | NA |  AU200 6.4 mm |
| B Series | Composi-Tight® Original | The original bands have a slightly shorter length. | | | | |
| | thickness 0.0015" |  B050 3.2 mm w/ext 4.2 mm |  B100 4.6 mm |  B150 3.8 mm w/ext 5.4 mm | NA |  B200 6.4 mm |
| MC Series | Composi-Tight® 3D Clear | Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands. | | | | |
| | thickness 0.0020" | NA |  MC100 5.0 mm | NA |  MC175 5.9 mm |  MC200 6.8 mm |

Item
FXH100
FXH150
FXH175
FXH200
FXH300

Item
FX100
FX150
FX175
FX200
FX300

Item
SXR100
SXR150
SXR175
SXR200
SXR300

Item
SM050
SM100
SM150
SM175
SM200
SM300

Item
M050
M100
M150
M175
M200
M300

Item
AU050
AU100
AU150
AU200
AU300

Item
B050
B100
B150
B200
B300

Item
MC100
MC175
MC200