

# Composi-Tight® 3D Fusion™

## Full Curve Non-Stick Matrices



**Wrap it up** - Full Curve bands wrap further around the tooth, out of the way, making ring placement easier. Still dead-soft and 0.0016" thin, they're the perfect choice.


**Keep tabs on 'em** - An occlusal placement Grab-Tab™ gives you excellent control with any placement instrument you choose. Just bend it over the adjacent tooth after you've positioned the band.



**Go deep or wide** - Improved subgingival extensions and added wrap-around length mean improved results on those more challenging restorations.

**Easy in** - Even though Composi-Tight® 3D Fusion™ Full Curve Matrices are made from the same ultra-thin dead-soft stainless steel as our other non-stick bands, their advanced curved form improves structural rigidity, making them our easiest-to-place band ever.

**Easy out** - Composi-Tight® 3D Fusion™ Full Curve matrices wouldn't be complete without the Slick Bands™ coating invented by Garrison. It reduces the adhesion of bonding agents and composites by a massive 92%. Regardless of your materials, these bands will slide right out!



### Composi-Tight® 3D Fusion™ Bands

**Description**

Assorted Band Kit: 5 sizes  
Assorted Band Kit: 5 sizes

**Item**

FXB01  
FXB02

**Refills**

Bicuspids  
Molars  
Bicuspids w/deep cervical preps  
Bicuspids and small molars  
Molars w/deep cervical preps

**Item**

FX100  
FX150  
FX175  
FX200  
FX300

\*Canadian MSRP

Be sure to check out the website at [www.garrisdental.com](http://www.garrisdental.com) for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.


Distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200



# Garrison BAND COMPARISON CHART

Dental Solutions



FX Series  
 FXH Series  
 SXR Series  
 B Series  
 M - SM - AU Series

	Pediatric	Bicuspid	Bicuspid with deep cervical preps	Bicuspid and small molars	Molars	Molars with deep cervical preps
FXH Series	Composi-Tight® 3D Fusion™ Firm matrix bands with eZ-Place™ technology for easier placement.					
	thickness 0.0015"	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm
FX Series	Composi-Tight® 3D Fusion™ Full Curve Bands have non-stick coating and an occlusal placement tab.					
	thickness 0.0016"	NA	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm
SXR Series	Slick Bands™ XR Matrices have non-stick coating and Grab-Tabs™ for easier placement.					
	thickness 0.0016"	NA	SXR100 4.6 mm	SXR150 3.8 mm w/ext 5.4 mm	SXR175 5.4 mm	SXR200 6.4 mm
SM Series	Slick Bands™ Matrices have non-stick coating.					
	thickness 0.0016"	SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm
M Series	Composi-Tight® bands are contoured in 3 dimensions to provide an accurate contact at the tooth's natural height of contour.					
	thickness 0.0015"	M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm
AU Series	The Gold bands are the same as the Composi-Tight® bands.					
	thickness 0.0015"	AU050 3.2 mm w/ext 4.2 mm	AU100 4.6 mm	AU150 3.8 mm w/ext 5.4 mm	NA	AU200 6.4 mm
B Series	The original bands have a slightly shorter length.					
	thickness 0.0015"	B050 3.2 mm w/ext 4.2 mm	B100 4.6 mm	B150 3.8 mm w/ext 5.4 mm	NA	B200 6.4 mm
MC Series	Clear bands are transparent and translucent to allow for trans-enamel polymerization without interference from metal bands.					
	thickness 0.0020"	NA	MC100 5.0 mm	NA	MC175 5.9 mm	MC200 6.8 mm

**Item**

FXH100  
 FXH150  
 FXH175  
 FXH200  
 FXH300

**Item**

FX100  
 FX150  
 FX175  
 FX200  
 FX300

**Item**

SXR100  
 SXR150  
 SXR175  
 SXR200  
 SXR300

**Item**

SM050  
 SM100  
 SM150  
 SM175  
 SM200  
 SM300

**Item**

M050  
 M100  
 M150  
 M175  
 M200  
 M300

**Item**

AU050  
 AU100  
 AU150  
 AU200  
 AU300

**Item**

B050  
 B100  
 B150  
 B200  
 B300

**Item**

MC100  
 MC175  
 MC200