

Rally™

Mini Polishers for Composite

Rally gets you to the finish line faster!

The composite finish line that is. With zero reprocessing and sterilization hassles and a lustrous diamond final gloss, Rally mini-polishers get you across the finish line and on to your next patient.

Get to the point - Every time!

Single patient Rally means that you get a clean, fresh, properly shaped polisher every time.

Mini, but mighty!

Rally's small size and tapered shaft create a truly agile polisher capable of reaching the smallest anatomical features without blocking your visibility.

A clean get-a-way!

You don't want to leave a trail of crumbs or spattered polishing paste on the track. Rally's diamond impregnated high shine final polishers won't crumble during use and don't require paste to create a beautiful, durable shine.


Give your mechanic the day off! (Maybe your assistant too!)

Rally is fully assembled - no little mandrels to lose track of, no discs to snap on, there's just nothing you need to mess with - including cleaning and sterilization.




Rally™

Mini Polishers for Composite




(Not actual size.)




Remove - Choose the Blue

Polishers for fast removal
and contouring



Smooth - Choose the Red

Polishers for smoothing
and silky shine



Shine - Choose the Gray

Polishers for ultra
high shine

Item	Description	Qty.
FPCK01	Composite Polisher Complete Kit	60 (10 each)

Item	Refills	
FPCC010	Blue coarse grit point composite polisher	30
FPCC040	Blue coarse grit cup composite polisher	30
FPCM010	Red medium grit point composite polisher	30
FPCM040	Red medium grit cup composite polisher	30
FPCF010	Gray fine grit point composite polisher	15
FPCF040	Gray fine grit cup composite polisher	15



Slow speed right angle latch shank. Use with water.



Do not exceed
maximum speed



Wear breathing
protection



Wear eye
protection



Use waterspray
minimum 50ml/min

*Canadian MSRP



Be sure to check out the website at www.garrisondental.com for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.

Distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200

CURION
dentistry, curated.