

# Composi-Tight® 3D Fusion™

## Ultra Adaptive Wedges



**Combining Soft-Face™ adaptive materials,**  
a firm inner core and advanced mechanical features to produce  
a wedge that truly works. Every time.

**Get a grip** The universal gripping block gives you absolute control of the wedging process whether you're using a cotton pliers, hemostat or "pin" tweezers. Use what you like - Composi-Tight® 3D Fusion™ has you covered.

**A real softie** Composi-Tight® 3D Fusion's Soft-Face™ over-mold allows the wedge to do what no other wedge can truly do - actually adapt to interproximal irregularities.

**Fins with serious function**

Soft retentive fins smoothly fold down during wedge insertion and then spring back when clear of the interproximal space. Wedge back-out is a thing of the past.

**Tough on the inside**

The firm inner skeleton allows for easy insertion and tooth separation just like a traditional wedge.

**No pain - all gain**

The upturned ski-like tip and soft radius edges glide over the rubber dam or papilla for trauma free insertion.

**Keep it clean (and dry!)**

Composi-Tight® 3D Fusion's patented design easily follows the contours of the teeth sealing things up nice and tight while preventing overhangs.



#### Description

Composi-Tight® 3D Fusion™ Wedge Kit  
Composi-Tight® 3D Fusion™ Wedge Kit

#### Item

FXK4  
FXK4-M

#### Refills

X-Small Yellow  
Small Blue  
Medium Orange  
Large Green

#### Item

FXYL-M  
FXBL-M  
FXOR-M  
FXGR-M

#### Item

FXYL  
FXBL  
FXOR  
FXGR



\*Canadian MSRP

Composi-Tight® 3D Fusion™ Wedges fuse two materials together to produce the perfect combination of adaptability and retention.

**Garrison**  
Dental Solutions

# Composi-Tight® 3D Fusion™ Ultra Adaptive Wedges

Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



## Composi-Tight® 3D Fusion™ Wedge Kits

| Description                         | Item   |      |
|-------------------------------------|--------|------|
| Composi-Tight® 3D Fusion™ Wedge Kit | FXK4   |      |
| Composi-Tight® 3D Fusion™ Wedge Kit | FXK4-M |      |
| Refills                             | Item   | Item |
| X-Small Yellow                      | FXYL-M | FXYL |
| Small Blue                          | FXBL-M | FXBL |
| Medium Orange                       | FXOR-M | FXOR |
| Large Green                         | FXGR-M | FXGR |



## WedgeWands®

WedgeWands® are designed to twist off the wand, but many prefer to leave it attached to help retract the cheek and for easy removal after the procedure is complete.

## WedgeWands®

| Description    | Item |      |
|----------------|------|------|
| WedgeWands Kit | WK4  |      |
| Refills        | Item | Item |
| X-Small Yellow | SWYL | WYL  |
| Small Blue     | SWBL | WBL  |
| Medium Orange  | SWOR | WOR  |
| Large Green    | SWGR | WGR  |


## A+Wedge® Controls Gingival Bleeding

- Controls gingival bleeding with an aluminum sulfate coating that has an astringent effect on the gingiva
- Improves clarity of restoration area
- Eliminates need to apply an astringent separately



## A+Wedge® Kits

| Description    | Item  |      |
|----------------|-------|------|
| A+Wedge Kit    | GWAK4 |      |
| Refills        | Item  | Item |
| X-Small Yellow | GWYL  | WYL  |
| Small Blue     | GWBL  | WBL  |
| Medium Orange  | GWAOR | WOR  |
| Large Green    | GWAGR | WGR  |



## G-Wedge®

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.

## G-Wedge™

| Description    | Item  |      |
|----------------|-------|------|
| G-Wedge Kit    | GWK4  |      |
| Refills        | Item  | Item |
| X-Small Yellow | GSWYL | GWYL |
| Small Blue     | GSWBL | GWBL |
| Medium Orange  | GSWOR | GWOR |
| Large Green    | GSWGR | GWGR |

## SoftWedge®

Soft Wedges™ compress during insertion, then swell to help fill voids to ensure that your preparations stay clean and dry—an absolute must for today's bonded restorations!



## Soft Wedge™

| Description    | Item |      |
|----------------|------|------|
| Soft Wedge Kit | WDK4 |      |
| Refills        | Item | Item |
| X-Small Yellow | WDYL | WYL  |
| Small Blue     | WDBL | WBL  |
| Medium Orange  | WDOR | WOR  |
| Large Green    | WDGR | WGR  |

\*Canadian MSRP

Be sure to check out the website at [www.garrisdental.com](http://www.garrisdental.com) for info on other innovative Garrison products, demo videos and our 6-month money-back guarantee.

Distributed by Curion

curion.ca | EN: 1.800.667.8811 | FR: 1.800.211.1200

