

# Product Description

With a calibre of experience in producing dental dams for leading dentistry institutions, Four D Rubber has developed the most advanced dental dam on the market; Isodam<sup>®</sup>

The primary function of Isodam<sup>®</sup> is to isolate specific areas of the mouth during procedures, reduce the risk of infection, and provide comfort to the patient.

The Isodam<sup>®</sup> provides exceptional stretch, strength, and tear-resistance while eliminating the potential for allergic reactions associated with latex products.

# **Key Features**

### **Superior Elasticity**

The polyisoprene material used to create the Isodam<sup>®</sup> allows for an incredible up to 1200% elongation that virtually eliminates tearing even with a punched dam.

## Size and Gauge Options

The Isodam<sup>®</sup> is available in medium and heavy gauges and 5" x 5" or 6" x 6" squares, ensuring the best non-latex dental dam option is available for any procedure.

## **Economy Package**

For those that use Isodam<sup>®</sup> on a regular basis, our economy package holds 5 times more dams than our standard package, creating less waste and greater savings.

Visit **www.fourdrubber.com** to learn more or request a sample



FDA registered, CE marked, ISO 9001. 2015 and ISO 13485. 2016 accredited North America Inquiries: naorders@fourdrubber.com All Other Inquiries: info@fourdrubber.com For professional use only Keep in a cool dry place Product of United Kingdon www.fourdrubber.com