

INSTRUMENT SELECTION DECISION TREE



*Created by Kathleen Bokrossy, RDH,
President of rdhu

CLIENT

Periodontitis

Periodontal Health

Stage I

Stage II

Stage III

Stage IV

Universal Kit

Area-Specific Kit

- Mirror
- Flexplorer
- Probe

- Mirror
- Flexplorer
- Probe
- Furcation Probe

- Mirror
- Flexplorer
- Probe

- SharpJack
- 4R/4L
- Anterior-Syntette*
- Syntette*

- M23
- Furcator
- GR 1/2*
- GR 15/16*
- GR 17/18*

- SharpJack
- 4R/4L
- BH 5/6

- SharpJack
- Anterior-Syntette*
- Syntette*

*Standard or mini

*Mini

*Standard or mini

DIAGNOSIS
TREATMENT

Streamline Essentials Series



Tel. 1.800.667.8811
www.biscocanada.com

Effectively access more areas with the least number of instruments – for each stage of periodontal involvement.

The Streamline Essentials Series offer simplified tray set-ups that maximize procedural efficiency and effectiveness. Each instrument kit in this series contains LM's winning combination of advanced sharpen-free tips* and the industry's highest-rated ergonomic handles.

*Exceptions include Diagnostic Instrument Kits and Furcator K/S

FEATURES AND BENEFITS:

- 100% sharpen-free tips
 - Unparalleled edge retention & superior sharpness
 - More thorough calculus removal with less hand strain
 - More comfortable scaling for both patient and clinician
 - Improved productivity & clinical results
 - Excellent tactile sensitivity
 - Exceptional control even in wet environments
 - Reduced waste of sterilization pouches and chemical indicators*
- *IPAC: Clinicians may no longer sharpen chairside. Must sharpen sterile instruments and then re-sterilize prior to use.

SUGGESTED # OF KITS	KIT TYPE	KIT CONTENTS	SKU #
8	DIAGNOSTIC	ServoE Cassette Mirror Handle Williams Probe (Stainless Steel) Flexplorer (Stainless Steel)	6552E 28ES (blue/red) 51ES 11F-12FES
6 - 8	AREA-SPECIFIC	ServoE Cassette SharpJack Syntette (Standard or Mini) Syntette Anterior (Standard or Mini)	6552E 315-316SDE 215-216SDE 275-276SDE
	UNIVERSAL	ServoE Cassette SharpJack Columbia 4R/4L Barnhart 5/6	6552E 315-316SDE 219-220SDE 227-228SDE
3	PERIO DIAGNOSTIC	ServoE Cassette Mirror Handle Flexplorer (Stainless Steel) Williams Probe (Stainless Steel) Furcation Probe (Stainless Steel)	6552E 28ES (blue/red) 11F-12FES 51ES 20B-21BXS
4	STAGE I	ServoE Cassette SharpJack Syntette (Standard or Mini) Syntette Anterior (Standard or Mini) Columbia 4R/4L	6552E 315-316SDE 215-216MSD 275-276MSD 219-220SDE
3	STAGE II-IV	ServoE Cassette M23 Furcator KS (Stainless Steel) Mini Gracey 1/2 Mini Gracey 15/16 Mini Gracey 17/18	6552E 313-314SDE 259-260ES 201-202MSD 245-246MSD 247-248MSD