

LMDarkDiamond™

LM

feel the
difference



PLEASANT AND PRECISE
SCULPTING OF COMPOSITE



Diamond-like feel & precision

Restorative instruments with LM Dark Diamond™ non-stick coating are ideal for placing and contouring composite resins. The diamond-like coating reduces stickiness and allows pleasant and precise restorative instrumentation.

- Exclusive diamond-like coating for non-stick placement and sculpting of composite resin
- Not sensitive to scratches due to hard surface
- Smooth surface for easy cleaning
- Dark-colored coating maximizes contrast against fillings and minimizes reflections
- Ergonomic LM-ErgoSense® technology delivers excellent precision and feel



Scan the QR codes for videos and more information

Instrument range



Burnisher

LMErgoSense®
LM 31-32DD ES

- LM Dark Diamond™ non-stick coating
- For packing and modeling composite



Burnisher-Plugger, ball-ended, 1.5–2.5 mm

LMErgoSense®
LM 471-473DD ES

- LM Dark Diamond™ non-stick coating
- For packing and modeling composite



Modeller MH

LMErgoSense®

LM 708-709DD ES

- LM Dark Diamond™ non-stick coating
- For packing large proximal fillings
- Cylindrical working end will not get stuck in composite during light-curing



Modeller MH Mini

LMErgoSense®

LM 707-708DD ES

- LM Dark Diamond™ non-stick coating
- For packing small and medium proximal fillings
- Cylindrical working end will not get stuck in composite during light-curing



Packer-Modeller IO

LMErgoSense®

LM 482-483DD ES

- LM Dark Diamond™ non-stick coating
- For packing and modeling composites
- Conical working end optimal for shaping occlusal surfaces



Plastic instrument, anterior

LMErgoSense®

LM 43-47DD ES

- LM Dark Diamond™ non-stick coating
- For placing, packing, and contouring composite



Plastic instrument, posterior

LMErgoSense®

LM 48-702DD ES

- LM Dark Diamond™ non-stick coating
- For placing, packing, and contouring composite



Mini Spatula

LMErgoSense®

LM 444-445DD ES

- LM Dark Diamond™ non-stick coating
- For placing and contouring composite

LM Dark Diamond™ instruments can be cleaned and sterilized like any other LM™ instrument.

LM Dental™

 LM-Instruments Oy

Norrbyn rantatie 8 | FI-21600 Parainen | Finland

Tel. +358 2 4546 400

info@lm-dental.com | www.lm-dental.com

Exclusively distributed by Curion



2571 Smith Street Richmond, BC., V6X 2J1 Canada

EN: 1.800.667.8811 | FR: 1.800.211.1200

curion.ca | wecare@curion.ca

