

LM

feel the difference

LM ErgoMix

Gentle titanium for implants **NEW!**

A great new instrument series for implant maintenance. Tips are made of sensitive but still effective titanium, which will not harm the abutments. The ErgoMix structure ensures excellent tactile sensitivity and optimal rigidity.



LM ErgoMix spares

Spare tips titanium	
Implant Mini Universal Curette L	LM 283MTi EM
Implant Mini Universal Curette R	LM 284MTi EM
Implant Mini Gracey 1	LM 201MTi EM
Implant Mini Gracey 2	LM 202MTi EM
Implant Mini Gracey 11	LM 211MTi EM
Implant Mini Gracey 12	LM 212MTi EM
Implant Mini Gracey 13	LM 213MTi EM
Implant Mini Gracey 14	LM 214MTi EM

Exclusively distributed by Curion

2571 Smith Street Richmond, BC., V6X 2J1 Canada
EN: 1.800.667.8811 | FR: 1.800.211.1200
curion.ca | wecare@curion.ca



LM

feel the difference

LM ErgoMix



Stay sharp & go green with exchangeable tips

Renowned LM ergonomics

Take the green way to sharpness

The new hand instrument series LM-ErgoMix offers the dental team a specially economical, ecological and ergonomic alternative. You can forget sharpening, but still work always with the best cutting edge. When ever you feel a need for a tip replacement, simply open the lock grip, insert a new tip and close the lock grip. No extra tools needed! Special alloyed LM-DuraGradeMAX super steel guarantees long lasting sharpness.



Lock grip – sealed and secure construction

Use LM-ErgoMix in the normal way, clean and sterilize it like any other instrument. The sealed construction guarantees high level of hygiene.



Retipping offers many benefits:



- always sharp instruments
- saves time from sharpening
- minimizes waste
- enables personalized tip combinations