


Now faster and more predictable!



Reliable, tight contacts
Drawn-wire nickel titanium plus PEEK produces consistent separating pressure longer!

Easier, more predictable contours
Marginal ridge enhancements built into the matrix bands AND rings help provide that finishing flair!

Faster, less frustrating Class IIs
The Strata-G™ rings, bands and wedges combine for a 270° Seal™ for ultimate flash control. The rings won't spring off the tooth, and the amazing Wide Preparation ring makes tough wide preps no problem.

Strata-G™ forms a 270° Seal™ for ultimate flash control.

the all new
STRATA-G™
order today!



- SGR-KSH-10** contains:
- 1 Short Ring (blue)
 - 1 Tall Ring (orange)
 - 1 Wide Prep Ring (green)
 - 50 Assorted Matrix Bands (10 each SG050, SG100, SG175, SG200, SG300)
 - 50 Assorted Firm Matrix Bands (10 each FXH100, FXH150, FXH175, FXH200, FXH300)
 - 10 VariStrip™
 - 80 Assorted Wedges (20 each SGYL, SGBL, SGOR, SGGR)
 - 12 Assorted Rally™ Polishers (2 each FPCC010, FPCC040, FPCM010, FPCM040, FPCF010, FPCF040)
 - 1 Ring Placement Forceps (FXP01)


Introductory Price
\$875 CAD MSRP



All Garrison products are covered by our six-month money-back guarantee. If you are not satisfied with your product, return the unused portion and we will issue a refund.

It's that simple!

the all new
STRATA-G™




Garrison.

Garrison® sectional matrix systems are the standard of care.


Tofflemire-style System

- Fails to restore proximal anatomy
- Thin contact at the marginal ridge
- Large food trap (below)
- Increased likelihood of fracture, occlusal interference, recurrent caries and periodontal disease





Strata-G™ Sectional Matrix System

- + Operator-friendly retaining system
- + Naturally contoured bands
- + Anatomically correct contacts
- + Contacts at the height of contour
- + Rings produce optimal tooth separation for excellent, tight contacts



270° Seal™ for ultimate flash control!



Longer lasting by design.

Advanced materials and design combine to produce our longest lasting, most predictable separator ring ever!

Strong tooth separation yet easier to open!


Nitinol drawn-wire nickel titanium is aligned at the molecular level for superior performance yet requires less effort to open. Garrison's proprietary StrataSet™ treatment maximizes the longevity of this advanced material.



More uses!

The PEEK ultra-durable reinforcement on the back of the ring allows Strata-G™ to retain a higher percentage of initial strength for much longer than non-reinforced rings. Proprietary StrataBond™ silicone/PEEK tips withstand hundreds of applications.

Improved anatomy—Unbeatable grip




Enhanced marginal ridge guides and supports matrix bands.

StrataBond™ silicone for unsurpassed flash control.

Ultra-Grip™ retention extensions provide tenacious grip and are gentler to the papilla.

Flexibility for a better cervical seal!

Strata-G™ Wedges utilize Soft-Face™ materials and advanced design to promote a better fit.



Stack them up

The smooth tunnel design easily accommodates stacking of wedges. This feature is extremely helpful in cases of gingival recession or when simultaneous wedging from the buccal and lingual is desired.

No back-out

Soft retentive fins smoothly fold down during wedge insertion and spring back when clear of the interproximal space, locking Strata-G™ Wedges in place.


Gentle on the papilla

Smooth tunnel design and rounded edges glide over the papilla with less trauma.



A better seal by design

The inverted V shape sits lower and seals deeper interproximally while flexible side walls and Soft-Face™ TPE over-mold allow Strata-G™ Wedges to easily adapt to concavities and irregularities.



Anatomy top to bottom

Five sizes of super-thin stainless steel in both dead-soft and firm ensure you have the best matrix band to restore every case.



Optimized marginal ridge

Strata-G™ Bands provide shape and guidance for developing the marginal ridge of Class II restorations without blocking your access.

ProfileBooster™ subgingival apron


Achieving a natural emergence profile in challenging cervical areas is easy with Strata-G™ Bands—it's built right in!



- Non-stick and color coded
- Easy placement tabs

- Dead-soft and only 0.0016" thin
- NEW pedodontic size!

Firm Bands for challenging cases



A critical addition in the Strata-G™ kit, Firm Bands are specially tempered stainless steel and only 0.0015" thin. Their stiffer nature makes them ideal for unbroken contacts and extra deep or extra wide preparations.

Wider indications—wider smiles!

Dentistry's premier Wide Preparation separator ring (green) makes Strata-G™ the most versatile and user-friendly system available. Tackle the most challenging posterior restorations with confidence knowing that you will achieve excellent contact and contour.



faster—more predictable—flash control
STRATA-G™

PMCA223SG