

Discover why Dentists are
making the switch
to the NEW VacuLUX mouthpieces:

VacuLUX

HVE ISOLATION SYSTEM

VacuLUX

ISOLATION MOUTHPIECES

Dryshield Zyris Isolite

 Universal Fit Designed to fit VacuLUX, Zyris & Dryshield systems	X	✓	X
 Autoclavable Save time & money by autoclaving with less prep	✓	✓ <small>Complicated Disassembly & Prep</small>	X <small>Simple-Sterile 1 Piece, Easy-Access Design</small>
 Comfortable Soft, protective overlap for patient's tongue & palate	X	✓	✓
 Mouthpiece Material Made of proprietary medical grade, optically clear silicone	X	✓	X
 Light Transmission Integrated light guide for targeted intraoral illumination	X	✓ <small>Optically Clear Engineered Silicone</small>	✓
 Quiet Calm Technology Engineered with open air pathways for 40% quieter procedures	X	✓	X
 Consistent Suction Evacuation chamber provides a continuous, uninterrupted dry field	✓	✓	X
 Liquid Rollover Edge Prevents pooling fluids for increased suction & spatter reduction	X	✓	X

**Simple &
Appordable**
Aerosol Control
& Isolation



Exclusively distributed by Curion
EN: 1.800.667.8811 | FR: 1.800.211.1200
curion.ca

Ascentcare
Dental Products

A system to fit your individual isolation requirements:



VacuLUX Mouthpieces

SKU: AVLUXA, AVLUXP, AVLUXS, AVLUXM, AVLUXL

Save time, money & the environment with these universal, autoclavable mouthpieces that will fit your current Dryshield, Zyris or VacuLUX system.

Designed with Simple-Sterile technology, limited nooks and easy-access cleaning flaps make sterilization fast, simple and safe.

Make the switch to the VacuLUX mouthpiece and your patients will notice a more comfortable and 40% quieter experience.

Made of proprietary, optically clear silicone for a longer lasting mouthpiece with optimal light transmission.

Available in four sizes: pedo, small, medium and large. Assorted packs also available.



VacuLUX HVE Adapter

SKU: VLUXD, VLUXZ

Want to expand your isolation abilities to every dental chair? Plug this aluminum adapter into any standard valve for instant isolation using Zyris, Dryshield or VacuLUX Mouthpieces.

Lightweight but durable aluminum construction will withstand repeated autoclave cycles. The one-piece design makes for simple cleanup.

A simple, portable & affordable isolation and aerosol control solution for your practice.



HyperLUX Light Engine

SKU: HLSTM

Illuminate your isolation mouthpiece and allow yourself the advantages of improved visibility.

The rechargeable, portable, and cost-effective light engine illuminates Dryshield, Zyris and VacuLUX mouthpieces.

Its patented, magnetically coupled fiber optic light guide allows for simple, plug & play illumination at any chair.

Operates at three settings: high, medium and low.

All VacuLUX adapters are equipped to allow illumination, so upgrade your system now or at any time.



Whisper Lite Kit

SKU: WLK5

This portable kit reduces sound and offers a more lightweight feel for your patient.

The Whisper Light kit works with any standard HVE valve in your operatory.

Kit includes a 5 foot kink-free vacuum hose, a stainless steel valve adapter tube, instrument resting clip and an aluminum coupler for the VacuLUX HVE adapter or any of your HVE dental instruments.

MicroValve

SKU: MCV

This lightweight, aluminum microvalve gives you convenient control of your adapter without the cumbersome weight & size of your standard valve.

Small, but mighty, this valve doesn't restrict airflow with its high-volume, large-bore 11.11mm inside diameter & can be used with any HVE instrument.