

August 13, 2019

Philips Respironics
1001 Murry Ridge Lane
Murrysville, PA 15668

RE: Airline Travel with DreamStation devices

To whom it may concern,

Philips Respironics has designed and tested the identified devices in support of the Code of Federal Regulations 14 CFR 382 “*Nondiscrimination on the Basis of Disability in Air Travel; Final Rule*”. Testing was completed per section 21, Category M, RTCA DO-160G EMI/RFI requirements and IEC 60601-1-2 4th edition, Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests.

Philips Respironics devices that are suitable for airline use have airline travel operating instructions included in the User manual. Additionally, the following symbol is located on the bottom of the device:



Complies with
RTCA/DO-160G
section 21, category M

The DreamStation User Manuals specify that the device is suitable for use on airlines.

Part Number	Device Description
AUX400S15	DreamStation CPAP Pro, AU
AUX400H15	DreamStation CPAP Pro w/ Humid, AU
AUX400T15	DreamStation CPAP w/ Humid/HT, AU
AUX500S15	DreamStation Auto CPAP, AU
AUX500H15	DreamStation Auto CPAP w/ Humid, AU
AUX500T15	DreamStation CPAP w/ Humid/HT, AU
AUX700S15	DreamStation Auto BiPAP, AU
AUX700H15	DreamStation Auto BiPAP w/ Humid, AU
AUX700T15	DreamStation CPAP w/ Humid/HT, AU
AUX400T15C	DreamStation CPAP Pro w/Hum/HT/Cell, AUS
AUX500T15C	DreamStation Auto CPAP w/Hum/HT/Cell, AUS
AUX700T15C	DreamStation Auto BiPAP/Hum/HT/Cell,AUS
AUX400S15	DreamStation CPAP Pro, AU
AUX400H15	DreamStation CPAP Pro w/ Humid, AU
AUX400T15	DreamStation CPAP w/ Humid/HT, AU
AUX500S15	DreamStation Auto CPAP, AU
AUX500H15	DreamStation Auto CPAP w/ Humid, AU
AUX500T15	DreamStation CPAP w/ Humid/HT, AU
AUX700S15	DreamStation Auto BiPAP, AU
AUX700H15	DreamStation Auto BiPAP w/ Humid, AU

If you have any other questions regarding our products, please feel free to call the Philips Respironics Customer Services department 1-724-387-4000 or +49 8152 93060. You can also use the following address:

Philips Respironics
1001 Murry Ridge Lane
Murrysville, PA 15668

Sincerely



8/13/19

Daria Brown
Sr. Regulatory Affairs Manager, SRC
1740 Golden Mile Highway
Monroeville, PA 15146