




AirFit™ N30
Nasal cradle mask

Simplicity from the start

ResMed AirFit™ N30 is the CPAP mask that discreetly fits into your lifestyle to help you stay on therapy. Our lightest mask ever* features a soft, curved nasal cradle cushion and an adjustable QuickFit™ elastic slip-on headgear.

Awaken your best fit.




- 1 QuickFit elastic headgear with split strap** is designed to make AirFit N30 easy to fit. It aims to provide stability and adjustability with minimal facial contact.
- 2 Single side clip** allows you to make minor adjustments to the headgear, if needed.
- 3 Nasal cradle cushion** rests discreetly under the nose and is curved for comfort.
- 4 QuietAir™ vent technology** is designed to break up exhaled air so it's spread out lightly, minimising sleep disruptions.
- 5 Slim short tube** is designed to be small, light, and flexible, offering freedom of movement without compromising seal.

**14
DAY**

**MONEY BACK
GUARANTEE**

We're so sure that you will love our AirFit N30 nasal cradle mask, that we are offering a 14 day money back guarantee if you are not completely satisfied.

Terms and conditions apply.



46 minutes more sleep

A ResMed study conducted in 2019 showed that on average, CPAP users slept 46 minutes longer each day with AirFit N30 compared to their current mask, AirFit P10.[†]

SIMILAR OPTIONS



AirFit™ P10

Nasal pillows mask

AirFit P10 is designed to be as comfortable and unobtrusive as possible, and offers premium comfort and ease of use without compromising stability.



AirFit™ N30i

Nasal cradle mask

AirFit N30i features a top-of-the-head tube design and nasal cradle cushion that sits comfortably under the nose, so you can sleep in many positions. A great option for side sleepers.



Speak to your ResMed Authorised Dealer for more information or to purchase the AirFit N30.

"It can't get any easier than this."[‡]

- AirFit N30 study participant

* As of October 2019. The total weight of AirFit N30, which includes the headgear and excludes the packaging, is 45g.

† Statistically significant differences ($p < 0.05$) in median scores were found between the two masks, AirFit N30 and AirFit P10.

‡ ResMed guided external clinical study of 17 evaluable CPAP mask users from 12–25 February, 2019 in Sydney, AU who trialled ResMed AirFit N30 at home for 7 nights in place of their prescribed mask, ResMed AirFit P10. CPAP users rated subjective criteria for AirFit N30 and AirFit P10 on a Likert Scale of 0–10.

AirFit, QuickFit and QuietAir are trademarks and/or registered trademarks of the ResMed family of companies. For patent information, see ResMed.com/ip. Specifications may change without notice.

© 2019 ResMed Ltd. 10111641 2019-10

