


# EcoStar Range

The perfect match for patients and providers looking for performance and simplicity in CPAP therapy

T R E A T M E N T



- ▷ Small, compact and very lightweight design
- ▷ Universal Power supply
- ▷ Compatible with the GoodKnight® H<sub>2</sub>O heated humidifier
- ▷ Easy to carry
- ▷ Maximum comfort and almost anywhere in the world without specific adaptation

**EcoStar**, **EcoStar Info** & **EcoStar Auto** are the successors of the world-famous GoodKnight® 420 G, S & E. These devices adapt easily to the user's daily routine and guarantee high-end performance, comfort and ease of use.


# EcoStar Range

Designed to improve compliance of Obstructive Sleep Apneas (OSA) patients treated by Positive Airways Pressure (PAP), the EcoStar product range provides the correct pressure and reliable performance, day after day.

The EcoStar product range uses high end components and has a very low failure rate, making these devices extremely reliable and more than ever the best «workhorses» among the low-cost CPAP devices.

**EcoStar** (CPAP) is about simplicity without compromising performance featuring unmatched pressure stability, time ramp, automatic leaks compensation, automatic altitude barometric and pressure adjustment.

With its basic configuration **EcoStar** (CPAP) is a very economic and efficient PAP device.



**EcoStar**

**EcoStar Info**

**EcoStar Auto**



**EcoStar Info** (CPAP) & **EcoStar Auto** (Auto-CPAP) are equipped with state-of-the-art technologies and algorithms to optimize user comfort and compliance : comfort mode, resistance compensation for any type of mask and tube, intelligent ramp, hands-free start (Auto-on) and Auto-off are standard.

**EcoStar Info** & **EcoStar Auto**'s algorithm detect inspiratory flow limitations, snoring and residual apneas and hypopneas. It differentiates between obstructive and central events and allows the device to deliver the optimum pressure at any time.

**EcoStar Info** & **EcoStar Auto** save the efficacy data on leaks and usage for one year in the device's memory, along with the real patient flow over the last 16 hours of use. These data are transferred to the DreamStar Analyze software and presented in the form of summary reports and graphs for simple and complete assessment of the treatment efficacy.

**EcoStar Info** & **EcoStar Auto** work with our OxyLink oximetry and connect to our SEFAM Connect telemonitoring platform.



For more information, please contact :

[www.sefam-medical.com](http://www.sefam-medical.com)



## Customer Service

T. +33 (0)3 83 44 77 40

F. +33 (0)3 83 44 99 87

[customerservice@sefam-medical.com](mailto:customerservice@sefam-medical.com)



**Manufacturer** : SEFAM - 144 avenue Charles de Gaulle - 92200 Neuilly-sur-Seine - FRANCE  
**Factory** : SEFAM - 10 allée Pelletier Doisy - 54600 Villers-Lès-Nancy - FRANCE



BR0149-03 / 10-2017

**Sefam**