

Heel & Ankle Protector




- Product provides immediate moisture control as well as cushioning support to heel area
- Extra layers of fabric allows for heel suspension
- Product will keep the skin at a regulated temperature which will eliminate the likelihood of skin breakdown
- Product will eliminate odors generally caused by excess humidity
- Dri-release fibers found inside the product will effectively attract and expel moisture 4 times faster than regular fabric.
- The 3-dimensional structure inside the product allows for air circulation inside the product
- Soft nylon exterior increases patient mobility
- Adjustable straps are soft against the skin
- Fabric Content: 90% Polyester 10% Cotton
- Washing temperature 160 – 167°F (71 – 75 °C)
- * Drying Temperature 220°F



SIZE CHART

PRODUCT	CALF CIRCUMFERENCE
<input type="checkbox"/> SMALL-MEDIUM (VHAPB-SM)	8"-15"
<input type="checkbox"/> MEDIUM-LARGE (VHAPB-ML)	16"+



TexMedico

www.texmedico.com

T: 1-506-855-9959 • F: 1-506-382-0079 • Toll Free: 1.877.448.4432 (CANADA/USA)
Dieppe, NB, Canada