


MOISTURE CONTROL

A Closer Look at the 3 Dimentional Structure

Inside the Ventopedic Heel & Ankle Offloading Boot

> 4-Way Cushioning support



Study results confirm minimal moisture absorption levels under warm and humid conditions.

	Suggested Fitting Guidelines				
	Size	Code	Calf Circumference	Shoe Size (Women)	Shoe Size (Men)
	Small	VHAOLB-S	8" - 12"	Up to 6	Up to 6
	Medium	VHAOLB-M	13" - 16"	7 to 9	7 to 9
	Large	VHAOLB-L	16" - 20"	10 +	10 +
	Bariatric	VHAOLB-B	20"+		

VENTOPEDIC PRODUCT CATEGORIES













Manufacture's Contact Information:



T: 1-506-855-9959 • F: 1-506-382-0079 Toll Free: 1.877.448.4432 (CANADA/USA)

Dieppe, NB, Canada www.texmedico.com

The effective moisture & pressure reduction solution.



Lateral Malleolus

4-way cushioning support under the heel and around the lateral and the medial malleolus removes pressure on these 3 critical pressure ulcer areas. Heel



Studies indicate that over **36% OF ALL PRESSURE ULCERS** develop on these 3 body parts

ADDITIONAL FEATURES

- The inner liner is soft against the skin and will reduce the possibility of shearing and/or friction
- · Controlled levels of hydrophobic and hydrophilic fibers found in the inner liner allows the fabric to remain dry even under extreme humidity
- 4-way cushioning support under calf, achilles heel and around foot area
- · Soft and breathable nylon exterior
- Machine washable: 167F (75C) Machine dryable: 220F
- Manufacture Guarantee of 12 months on material defects and workmanship under normal usage conditions for the 12 month period
- · ID tag for patient information
- Latex Free
- 3 adjustable straps are soft and durable even after repeat washings
- Available in Small, Medium, Large and Bariatric sizes