



## **Southwest Medical Refurbishing Process for Anesthesia Systems**

- **Run Serial Number With OEM For Necessary Upgrades, Recalls, And Revisions**
- **Verify Presence, Function & Performance Of Internal Components Including, But Not Limited To:**
  - o **Pressure Gauges**
  - o **Flow Head**
  - o **Processor & Program**
  - o **ORMc**
  - o **O2 Monitor**
  - o **Power Supply**
  - o **Battery Voltages**
  - o **Alarm Channel**
  - o **Ventilator**
- **Repair, Rebuild, Or Replace Any Internal Components Listed That Do Not Meet OEM Specifications**
- **Perform Complete System Disassembly If Necessary**
- **Disinfect ALL Components**
- **Sand & Repaint All Previously Painted Surfaces Including Frame Drawers & Flow Head**
- **Resurface Tabletop To Original Appearance By 5-Step Sanding & Polishing Process**
- **Replace CTRs If Necessary**
- **Clean & Rebuild Absorber**
- **Reassemble Refurbished Components**
- **Install New Latex-Free Hoses & Bellows If Applicable**
- **Perform Comprehensive Manufacturer PM Service Prior To Staging**
- **Configure System To Customer's Specified Parameters**
- **Stage System For Optimum Performance & Reliability**
- **Label Machine & Components To Specified Requirements**
- **Perform Comprehensive Manufacturer PM Service Upon Delivery To Customer**

**Call or Email Southwest Medical for Additional Information**

**(800) 259-5940**

**[www.swmedical.com](http://www.swmedical.com)**