



Southwest Medical Refurbishing Process for Stretchers

- ***Run Serial Number With OEM For Necessary Upgrades, Recalls, And Revisions***
- ***Disassemble And Disinfect All Components***
- ***Polish, Replate, Repaint, Or Replace All Exterior Components***
- ***Install Rebuild Kits In All Hydraulic Jacks***
- ***Inspect All Working Components For Wear & Replace Defective Parts As Necessary***
- ***Soak, Degrease, Clean & Lubricate Wheel & Caster Assembly***
- ***Replace Foot Pedals If Necessary***
- ***Sand, Prime & Paint Frame To Meet New Specifications***
- ***Restore Stretcher Base & Lift System To Original Mechanical Specifications***
- ***Recalibrate Brake & Steer System***
- ***Reassemble Refurbished Components***
- ***Inspect And Lubricate All Pivot Points***
- ***Inspect And Lubricate Instant CPR Cables And Latch Release (If Applicable)***
- ***Inspect And Lubricate Siderail Latches***
- ***Apply New Decals And Safety Information***
- ***Stage Test Weight Tolerances And Table Movement Through All Available Articulations***

Call or Email Southwest Medical for Additional Information

(800) 259-5940

www.swmedical.com