



## ***Southwest Medical Refurbishing Process for Surgical Tables***

- ***Run Serial Number with OEM for Necessary Upgrades, Recalls, and Revisions***
- ***Disassemble and Disinfect all Components***
- ***Polish, Replate, Repaint, or Replace all Exterior Components***
- ***Clean, Test, and Recalibrate all Motors***
- ***Inspect, Clean, Acid Dip, and Polish all Working Components for Reuse***
- ***Remove and Ship Electrohydraulic Motor to Outsource for Refurbishment***
- ***Replace Main lift Cylinder as Needed***
- ***Sand, Prime, and Paint Table Frame to Meet New Specifications***
- ***Restore Table Base and Pedestal to Original Specifications***
- ***Replace Casters and Wheels as Needed***
- ***Restore Table Gear Train with Combination of New and Refurbished Gears, Bearings, and Clutches***
- ***Replace all Bushing Eccentrics (if Applicable)***
- ***Resurface Stainless Steel to Original Appearance by 5-Step Sanding and Polishing Process***
- ***Reassemble Refurbished Components***
- ***Apply New Decals and Safety Information***
- ***Stage Test Weight Tolerances and Table Movement through All Available Articulations***

**Call or Email Southwest Medical for Additional Information**

**(800) 259-5940**

**[www.swmedical.com](http://www.swmedical.com)**