

Two Year Warranty from LEAP • No Bungee Cords • 2D Barcode Reader • Higher Capacity


RTC Robotic Tool Changer Sample Loader
Sample Prep Automation, Push-Button Operation

RSI Robotic Sample Injector Frontend Automation
High Capacity, High Productivity

LEAP's PAL Sample Control Software

- controls many setups, methods
- integrates with most chromatography & mass spec data systems
- high efficiency sample prep
- scheduling software is standard
 overlapped -looking-ahead sample prep and inject

or - Chronos Software



PAL3 for GC

- unattended method and syringe change
- liquid injection headspace extraction
- SPME
- 2D barcode reader
- cold needle inject
- liquid / liquid extractions
- sample prep

PAL3 for LC and Work Stations

- many sample prep steps possible
 - -dilution, reactions, filtration, reformatting
 - -from test tubes to 1536 well plates
- different syringe sizes and types
- high speed vortexer
- 2D barcode reader
- liquid / liquid extractions
- sample prep
- multi-valve capability
- ITSP sample prep
- low-to-no carryover





Gazelle18 UHPLC Pump/PAL-Bundle

Replace out-dated unreliable LC Front Ends
Optimizes I C-MS productivity, solvent bottle

Optimizes LC-MS productivity, solvent bottles to column heater. The pump has broadest range of pressure and flow rate in the industry. One driver controls two binary pumps. Chronos Software and Analyst™ Datasystem.