







## 解熱鎮痛剤の分離 2

### Analgesics 2

Sunniest C18 5 $\mu$ m, 250 x 4.6 mm i.d.



Column: Sunniest C18 5 $\mu$ m, 250 x 4.6 mm i.d.

Mobile phase:

A) CH<sub>3</sub>CN/20mM Phosphoric acid = 45/55

B) CH<sub>3</sub>CN/0.1% Formic acid = 45/55

Flow rate: 1.0 mL/min

Temperature: 40 °C

Detection: UV@230 nm

Injection volume: 10  $\mu$ L in B mobile phase