

We're changing
the face of sleep

PHILIPS

RESPIRONICS

sense **and** simplicity

Wisp

Performance has a new style



Developing solutions that meet the needs of a rapidly changing industry calls for fresh thinking and collaboration. At Philips Respironics we traveled around the world asking clinicians, homecare providers, and patients how we could improve their experience.



Clinicians wanted us to resolve the tradeoff between nasal and pillows masks, homecare providers wanted a solution that could reduce operational complexity by getting more patients served from one package, and patients wanted a better overall mask experience that improved their lifestyle.

When you put it all together, the answer is Wisp.

- Compact design comfortably fits more than **98%** of patients



- Patients said they had a **better quality of sleep** with Wisp over the leading pillows and nasal masks



Allies
in Better Sleep
and Breathing

At Philips Respironics we're listening carefully to you and working toward solutions that contribute to healthy patients and healthy businesses.



- Open field of vision for a feeling of **greater independence** — patients can read and even wear their glasses



Wisp delivers


During our CPAP patient preference and RT/ST ease of use studies, Wisp was preferred over the leading nasal and nasal pillows masks in these categories:

Patient preference study results	Comfort
	Ease of use
	Seal
	Quality of sleep
	Less nostril irritation

RT/ST preference study results	Fit & seal
	Ease of assembly & disassembly
	Ease of patient setup

- The minimized design and headgear clips make it **easier** for patients to use and clean

Resolving the tradeoff



Self-adjusting top strap
with tube holder loop

Reversible fabric frame —
soft suedette and silky
sateen — also available in
a clear silicone

Easy to adjust with simple
headgear clips



Tip-of-the-nose
cushion design

Elbow quick
release tabs

Our unique **tip-of-the-nose** cushion is designed to sit below the nose bridge for a minimal contact experience and superior seal.



Petite



Small/medium



Large

Each package contains all three cushion sizes, which fit on the same frame.

Order Wisp now to see how we're changing the face of sleep for you and your patients

Item	Part number	U.S. HCPCS code
Wisp, fabric frame		
Wisp mask with fabric frame, with headgear	1094051	A7034, A7035
Wisp mask with fabric frame, without headgear	1094056	A7034
Wisp, clear frame		
Wisp mask with clear frame, with headgear	1094050	A7034,A7035
Wisp mask with clear frame, without headgear	1094055	A7034
Replacement parts		
Petite cushion	1094086	A7032
Small/medium cushion	1094087	A7032
Large cushion	1094088	A7032
Clear frame	1094080	N/A
Fabric frame	1094081	N/A
Wisp headgear	1094082	A7035
Wisp headgear, large size	1094078	A7035
Wisp headgear clips	1094083	N/A

Manage your patient resupply with one solution, Fit for Life, included with every Wisp mask.

Pick your Wisp

Why patients preferred each version from an appeal perspective



Clear frame preference

- Silicone is less visible and doesn't draw attention to the face
- Translucent, looks smaller



Fabric frame preference

- Feels softer and more pliable
- Gets you closer to the feeling of your pillow
- Easier to sleep on your side

**Philips Healthcare is part of
Royal Philips Electronics**

How to reach us

www.philips.com/healthcare
healthcare@philips.com

Asia
+49 7031 463 2254

Europe, Middle East, Africa
+49 7031 463 2254

Latin America
+55 11 2125 0744

North America
+1 425 487 7000
800 285 5585 (toll free, US only)

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)

Philips Respironics International
Headquarters
+33 1 47 28 30 82

Philips Respironics Asia Pacific
+65 6882 5282

Philips Respironics Australia
+61 (2) 9947 0440
1300 766 488 (toll free, Australia only)

Philips Respironics China
+86 400 828 6665
+86 800 828 6665

Philips Respironics Deutschland
+49 8152 93 06 0

Philips Respironics France
+33 2 51 89 36 00

Philips Respironics Italy
+39 039 203 1

Philips Respironics Sweden
+46 8 120 45 900

Philips Respironics Switzerland
+41 6 27 45 17 50

Philips Respironics United Kingdom
+44 800 1300 845

www.philips.com/respironics

Respironics and Wisp are trademarks of Respironics, Inc. and its affiliates.
All rights reserved.

Please visit www.philips.us/wisp



© 2012 Koninklijke Philips Electronics N.V.
All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from use of this publication.

CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.
Hoech RW 9/20/12 MCI 4105239 PN 1103630