

FREAKYS COLFAX - #1
5701 E. Colfax
Denver, Colorado
303.765.2628

FREAKYS NORTH DENVER - #2
5326 W. 41st Avenue
Denver, Colorado
303.458.7326

FREAKYS THORNTON - #3
9205 N. Washington
Thornton, Colorado
303.227.1535

FREAKYS GLENDALE - #4
4521 Leetsdale Drive
Glendale, Colorado 80222
303.329.0420

FREAKYS PLATTE #5
308 E. Platte
Colorado Springs, Colorado
719.536.0400

FREAKYS ACADEMY - #6
1714 Brookwood
Colorado Springs, Colorado
719.575.0420

FREAKYS SANTE-FE - #7
776 Sante Fe Drive
Denver, Colorado
303.413.0420

FREAKYS LAKEWOOD - #8
3355 S Wadsworth H-101
Lakewood, Colorado
303.736.2130

FREAKYS N. LAKEWOOD - #9
155 N Wadsworth
Lakewood Colorado
303.237.0420

FREAKYS PARKER - #10
10471 S PARKER ROAD
Parker, Colorado
303-955-8737

FREAKYS ARVADA - #11
10050 RALSTON ROAD
Arvada, Colorado
720-420-9377

FREAKYS LONE TREE - #12
8850 Maximus Drive
Lone Tree, Colorado
720.379.5523

FREAKYS WAREHOUSE - #101
720.541.7451



FREAKY'S SMOKE SHOP PRODUCT TRAINING MANUAL & HANDBOOK

Table of Contents

NOTES:

PRODUCT TRAINING

DETOX	4
DEFINITIONS	16
BANGARS	19

GENERAL:

NOTES	48
-------	----

NOTES:

710 Premium Grade Detox

There's nothing stronger! 710 Detox is the most effective one-hour detox regime on the market today. Formulated to remove toxins from concentrates and edibles, this highly specialized product is made for today's modern lifestyle. This fast-acting, one hour liquid drinks contain a highly specialized combination of powerful cleansing herbs and essential vitamins that work with your body's natural detoxification system to remove harmful toxins. Regardless of the level of toxins in your system this product will assist your body in eliminating unwanted toxins.

For over 20 years millions of satisfied customers have enjoyed the peace of mind that comes with using our **Ready To Drink** formulas for people with immediate cleansing needs, normal to moderate toxin levels and weigh up to 240 pounds.

PURE-DETOX 710 Professional Grade Detox, Directions For Use:

1. Pick a day for using 710 Detox
2. Avoid exposure to toxins for at least 24 to 48 hours before using 710 Detox
3. Shake bottles well as some settling may have occurred.
4. Drink one bottle of 710 Detox and consume the entire packet of 6 capsules.
5. Immediately refill the empty bottle with water, shake well and drink
6. Wait 15 minutes and consume the second bottle of 710 Detox
7. Immediately refill the empty bottle with water, shake well and drink
8. Urinate frequently after drinking 710 Detox to expel toxins
9. 710 Detox is most effective after three to four normal urinations

Get the Most Out of Your Cleanse with These Important Tips:

1. Drink as much water as possible on a daily basis.
2. Avoid taking OTC drugs large quantities of vitamins and acidic liquids such as vinegar and juices.

NOTES:

NOTES:

3. Avoid nicotine and caffeine and other unwanted toxins prior to use of the product.
4. Do not eat large meals before using the product.
5. 710 Detox has been formulated as part of an ongoing detoxification and internal cleansing program that includes a healthy diet and regular exercise.
6. Drink 6 to 8 twelve ounce serving of water per day for at least two days prior to starting your detoxification program.
7. During your detoxification program avoid foods that are high in saturated fat.
8. Avoid soda pop, fruit juices and coffee.
9. Avoid alcohol, over-the-counter drugs and other medications when possible.
10. Eat light meals including fresh fruits, vegetables and foods with high fiber content.
11. Women who are pregnant or breast-feeding should consult with their physician before using 710 Detox.
12. 710 Detox is not recommended for children.

Once you've completed these easy steps simply relax...you have purchased the best cleansing product on the market. Our formulations and products have been used by millions of people who are serious about cleansing and detoxification.

PURE-DETOX Extra Strength

<http://www.puredtx.com/wp-content/uploads/2012/04/3D-20oz-Pineapple-Orange-group.jpg>

Product Description:

This potent easy to use liquid drink is formulated to cleanse and detoxify men and women just like you. For over 20 years millions of satisfied customers have enjoyed the peace of mind that comes with using our **Ready To Drink** formulas for people with immediate cleansing needs, moderate to high toxin levels and weigh up to 240 pounds.

PURE-DETOX Regular Strength Directions For Use:

Follow these steps to ensure your satisfaction:

- Drink the entire contents of the bottle. Refill with water and drink at a comfortable yet consistent pace.
- PURE-DETOX starts to work immediately. Your cleansing and detoxification process has begun!
- Be sure to urinate frequently. This is very important to a successful cleansing program.
- If you desire, the additional consumption of water will extend your cleansing and detoxification regime for several hours.

Get the Most Out of Your Cleanse with These Important Tips:

- Drink as much water as possible on a daily basis.
- Avoid taking OTC drugs, large quantities of vitamins and acidic liquids such as vinegar and juices.
- Avoid nicotine and caffeine and other unwanted toxins prior to use of the product.
- Do not eat large meals before using the product.
- Once you've completed these easy steps simply relax...you have purchased the best cleansing product on the market. Our formulations and products have been used by millions of people who are serious about cleansing and detoxification.

Pure Detox Regular:

Follow these steps to ensure your satisfaction:

- Drink the entire contents of the bottle. Refill with water and drink at a comfortable yet consistent pace.
- PURE-DETOX starts to work immediately. Your cleansing and detoxification process has begun!
- Be sure to urinate frequently. This is very important to a successful cleansing program.
- If you desire, the additional consumption of water will extend your cleansing and detoxification regime for several hours.



Thermal Banger –

The fixed dish in the center helps retain heat, great for flavorful low temp dabs. Normally comes with a bubble cap.



Terp Slurper/ Vacuum Banger –

- It allows you to **use** large dabs or globs. ...
- It is perfect for low temperature dabbing as the concentrate comes into contact with the hottest part of the nail automatically without the need for a directional carb cap.

A **terp slurper** is a new design on an old method of dabbing. While buckets have been common for several years, this style of bucket helps ensure less wasted product regardless of the size of the dab consumed.



Recycler Banger –

Recycling power is no longer just for the dab rig itself, now your actual banger can start the recycling process right from the start. This feature helps make sure that you can use up every last bit of your wax, which will heat up more evenly as it circulates through this banger's recycling system.



Terp Spinner –

The thick bottom retains heat so you can take low temp dabs, the terp pearls, help evenly distribute your wax on the bottom of the banger while the cap spins the terp pearls.



Cadmium Banger –

this version contains a quartz core with only a cube of color changing material to indicate nail temperature. The colored cube inside the core will change to orange and red to let you know when the nail is the perfect temperature. The flat top design allows you to use a variety of different carb caps. This dab nail is a must have for any true dabbing enthusiast.



Diamond Banger –

With low temperature dabs being king, the fact that this design allows heat to travel and spread more evenly means that lower temperatures can be maintained in the **diamond** knot longer. It also ensures that the inside of the **banger** is a consistent temperature instead of having a few hot spots

Get the Most Out of Your Cleanse with These Important Tips:

- Drink as much water as possible on a daily basis.
 - Avoid taking OTC drugs large quantities of vitamins and acidic liquids such as vinegar and juices.
 - Avoid nicotine and caffeine and other unwanted toxins prior to use of the product.
 - Do not eat large meals before using the product
- One hour detox drinks are a great way to eliminate unwanted toxins from your body. Whether you have a long day of hiking ahead of you, are planning to run a marathon, perhaps you have low energy levels or have been stressed out lately and are looking to refresh your system. The reasons to detoxify are endless... the solution is simple. A fast acting, safe, high quality Herbal Cleanser.
 - A good Herbal Cleanser should contain of a variety of herbs, vitamins and minerals that are scientifically formulated to assist the body in a quick and effective detox. Herbal extracts such as Dandelion Root, Milk Thistle and Burdock Root are critical cleansing the liver and other vital organs of unwanted toxins.

Effective elimination of toxins:

- Eliminating toxins is what detox is all about. Our modern lifestyle is full of elements that create toxic build up in our body. Lack of exercise, fast food, smoking, drinking, environmental pollution – topped off with high stress and not enough sleep it's easy to see how the body can become overwhelmed with toxins that may cause a host of health issues. In today's highly toxic world toxins may build up faster than the body can eliminate them, this can lead to illness and chronic fatigue. Neometrx products assist the body in detoxifying and eliminating toxic waste from the body quickly and effectively. It's truly amazing how much cleaner and healthier you feel when your body has been properly detoxified.

Not all detox products are the same:

- Always choose a one hour detox that is a true herbal cleansing product, these can cost a little more but are well worth the investment. The cheaply made and less expensive one hour drinks are usually nothing more than fortified water containing only sugar, flavoring and vitamin B. The results can be rather disappointing.
- Make sure to check the labels and choose a product that will do the job right the first time.

PURE-DETOX Regular Strength Directions For Use:

Follow these steps to ensure your satisfaction:

- Drink the entire contents of the bottle. Refill with water and drink at a comfortable yet consistent pace.

- PURE-DETOX starts to work immediately. Your cleansing and detoxification process has begun!
- Be sure to urinate frequently. This is very important to a successful cleansing program.
- If you desire, the additional consumption of water will extend your cleansing and detoxification regime for several hours.

Get the Most Out of Your Cleanse with These Important Tips:

- Drink as much water as possible on a daily basis.
- Avoid taking OTC drugs large quantities of vitamins and acidic liquids such as vinegar and juices.
- Avoid nicotine and caffeine and other unwanted toxins prior to use of the product.
- Do not eat large meals before using the product.
- Once you've completed these easy steps simply relax...you have purchased the best cleansing product on the market. Our formulations and products have been used by millions of people who are serious about cleansing and detoxification.

PURE-DETOX Maximum Strength – Jamaican Punch:

- We packed this drink with all the nutrition, vitamins and herbs to cleanse and detox the body. We also added a melody of all-natural fruit flavors to really get the party going in your mouth. How awesome is it that you can detox your body with such a delicious tasting drink... You'll wish you were sitting on a beach while drinking it.
- http://www.puredtx.com/wp-content/uploads/2012/05/supp.fact_jam_punch_.png

HOW TO CHOOSE A DETOX FOR A CUSTOMER

- With so many detox products on the market today choosing the right detox solution for your specific needs can be a confusing and down right daunting task. Knowing a few key facts can help you to quickly and easily choose the right product for you.
- **Before you start looking for a product, ask yourself these questions: What are my detox goals? What part of my body do I want to cleanse? How soon would I like to see results?**
- Are you are looking for a fast acting-one hour detox that focuses on heavy metals and toxins ?
- Is your liver, kidney's and gallbladder your main concern?
- PURE-DETOX is formulated for goal oriented people, on the go, that are concerned about moderate to high levels of toxic build up in their bodies.
- **When you need something quick, fast acting liquids are the smart choice. Our powerful and effective cleansing drinks address 5 vital areas of detoxification that are critical for success.**

Bangers



Reactor Banger –

Ideal for low temp dabs, normally comes with a claw carb cap. The fixed thick piece of quartz in the center retains heat so you can get a better low temp dab.



Color Bottom Banger –

The large hollow nail on the bottom allows you to heat up a massive dose of oils and waxes. The thick colored base boosts heat retention and also gives the piece an elegant look.



Dome-less quartz banger –

domeless nail design that features a direct inject design, old school way of dabbing.

The technology behind the Health Stone relies on gravity and a bit of skill. It is made from a vapor stone with many tiny holes and gaps that allows the concentrate to sink into the stone and eventually vaporize.

First, you load your dab on the coin-sized, couple-centimeters-thick stone, keep your torch far away and melt your concentrate slightly into the stone. Then bring the flame closer, all while providing a strong consistent drag to vaporize the concentrate. Depending on the piece you're hitting, it could be difficult to have enough air to pull the vapor into your lungs.

It takes time to learn to hit a health stone correctly, and depending on the type you're hitting (aka whether or not you're dry hitting it) it could feel really harsh. This fact alone deters many people away from a health stone.

But if it's not this issue, it's the fact that you do inhale a percentage of the gas emitted from your torch. That being said, there are differing opinions on whether or not this should cause concern

The Nectar Collector

Last, but certainly not least in the clearly wide array of options is the [Nectar Collector](#). Now, this option does require the torch but provides a better portability factor than normal rigs do. The nectar collector is a vapor straw with the added benefit of water filtration.

Nectar collectors come with different types of tips, and are made from borosilicate glass. To use one you need a torch and a glass pad (preferably with a lip) to place your dab on.

What is alluring about this is the water filtration provided (which the handvapes that I'm familiar with don't have). Although this doesn't appear to necessarily make dabbing on the go any faster, it certainly looks like a way to vaporize on the go more efficient.

1. **Cleanses and detoxifies the liver, urinary tract and other vital organs.**
 2. **Stimulates the metabolism.**
 3. **Purifies the blood.**
 4. **Boosts the immune system.**
 5. **Fiber to bind and draw toxins and heavy metals out of the body.**
- **Try PURE-DETOX for yourself...Millions of people have depended on our formulations and products for over 20 years. When results count our formulas work fast, the first time.**

DETOXIFY GREEN CLEAN

http://www.detoxify.com/wp-content/uploads/2013/08/GreenClean_Box_39545_r05_80C_1000px-SQ.jpg

- Detoxify Green Clean is the latest innovation in detox drinks. The two 4oz bottles contain the most powerful blend of cleansing herbs and metabolic boosters on the market today. No detox drink is more powerful, better tasting, or easier to use – Honey Tea Flavor. No need to PreCleanse!
- You can enjoy a healthier, cleaner body with Detoxify Green Clean. Green Clean is the herbal cleansing detox drink concentrate used by people just like you. They cleansed their systems of toxins and impurities with Green Clean, and so can you.
- Detoxify Green Clean is the way to help lower the impurities your body absorbs every day. Simply follow these instructions on how to use Green Clean:

Detoxify Green Clean Directions:

- Pick your day for using your Green Clean for intensive cleansing.
- Begin your cleansing program with Green Clean.
- Shake the Green Clean Cleansing Blend well and drink entire contents of the bottle. Next, shake Green Clean Green Tea Metaboost well and drink entire contents of bottle.
- It is permissible to mix both Green Clean bottles with up to 16oz of water.
- Wait 15 minutes. Drink 32oz of water.
- Continue to drink plenty of water throughout the day to extend your cleansing program for hours.
- Urinate frequently. Frequent urination indicates that you are experiencing optimal cleansing with Green Clean.

Additional Information:

- Green Clean is intended for periodic intensive cleansing that is a part of an ongoing cleansing routine.
- For best results, use Detoxify's PreCleansing products for 2 days prior to using Green Clean.
- For optimal cleansing and health, take Detoxify Constant Cleanse every day with six 16oz glasses of water.
- Eat light meals, including fruits, vegetables, and fiber during your cleansing program

DETOXIFY XXTRA CLEAN

<http://www.detoxify.com/wp-content/uploads/2014/01/>

small layer of titanium oxide. If this layer is especially thick from daily use, it can effect the nail's productivity.

Of course, your nail will become seasoned over time so some people bypass this step, but the likelihood your dabs will taste like metal for a while is high. These nails come in domeless versions as well.

Quartz Nails

The next type of nail on the list is the Quartz nail. These nails are considered to reign supreme for savoring the [terpenes](#) in concentrates, meaning that you can get a really amazing taste.

Quartz nails are similar to glass in that they take a little longer to heat than the Ti nail. However, they are much more affordable in comparison.

Just like their glass cousin, quartz nails are prone to cracking from extensive use. But quartz nails can last you a very long time with care.

Ceramic Nails

What we do know about the ceramic nail is that it's non porous, so you don't need to season it. Also, it is said to be chemically inert, making it medical/ food grade product.

Many people who have tried the ceramic nails say that the taste you get is great, but it's nearly impossible not to break or crack from daily use. . Even the most careful users could still potentially see cracks or fractures.

Another downside about ceramic nails for the average dabbler is that the heating time takes considerably longer than Ti, glass or quartz. They require patients to heat the nail for approx. 2 minutes before and to let it cool down for about a minute before dabbing.

The upside to this heating process is that ceramic is remarkable for retaining heat, so it's possible to take more than one dab from heating it once – which could be a benefit to some.

Skillets, Curves, and Honey Buckets

Instead of the more traditional nail, some people prefer the use of a skillet. How that works is you heat a Ti pad attached to a swing arm to the desired temp, swing it under a bell or curve, and apply your dab while you inhale through your rig. This often results in vapor loss, and has been deemed unfavorable by many.

Similar to a curve is a honey bucket, or honey pot. Whats appealing about the concept of a honey bucket is that there is no exposed nail like there is with of any type domeless nail, but no ability to lose the dome for a nail that not only requires one, but is attached to it.

With the swing arm attachment, you eliminate a lot more risk of burning yourself or something else since the nail is protected by the rest of dome.

The Health Stone

What makes the [Health Stone](#) stand out is that it's the quickest way apart from vapor pens to dab. All you need is a triple flame torch and some oil – no need to wait for a nail to be heated, and no domes are required.

The Dabber

The first utensil needed is a tool to transport your concentrate to the nail – this is referred to as a dabber. Dabbers can be made out of anything from glass to titanium (Ti) to stainless steel. It's important to remember that it should not melt into the nail when you apply your dab. The dabber also should not affect the taste of your concentrate.

Interestingly, there are specific types of dabbers that work best with different consistencies of concentrates. It's been said that dabbers with a ballpoint work best with oily, stickier concentrates, while flat heads work best to cut and break into shatter.

If you know which type of concentrate works best for you, you will soon find the best dabber to use it with.

CARB CAPS

Carb caps basically give your dab a one-two punch during vaporization at a lower temp.

To explain, the carb cap combines both conduction-style and convection-style heating on a domeless nail by covering it with a cap. This allows your concentrate to spread out and be evenly vaporized, while having the added benefit of allowing you to dab at lower temperatures.

Not all nails are capable of using these carb caps, but some companies have made their own to fit their specific product since they emerged.

Glass, Ti, or quartz nails often require domes to keep the vapor from escaping, making aking it impossible to use a carb cap.

Glass Nails

Most people start off with the glass nail that commonly comes with a new concentrate rig, and they can be effective. However, glass nails are prone to heat cracking over time.

Personally, I've had my glass nail for a long time, and like to keep it in the rotation – I like the added fact that it's a natural alternative to Ti.

All in all, glass nails do the trick while in a pinch. However, it can be difficult to achieve a satisfactory dab each and every time.

Titanium Nails

The Titanium (Ti) nail is the average dabber's go-to method of choice. Most smoke shops will recommend the Ti nail over other nails due to the fact that it takes less time to heat compared to the glass. They are generally more reliable in terms of daily use, but they can often disrupt the flavor of your concentrate.

To best ensure the taste is as accurate as possible, you must season your nail. They recommend seasoning with oil, not any other method (some prefer to use water).

Although these nails are the most user-friendly in terms of quickness, they aren't the best for patients in comparison to glass, quartz, or ceramic competitors. For instance, your Ti nail will oxidize over time, causing a

XXClean Tropical Fruit Box-Bottle 1000pxh_sq-300x300.jpg

- Detoxify Xxtra Clean is the cleansing herbal detox drink formulated for people with higher toxin levels and time to PreCleanse. The 20oz ready-to-drink formula is packed with vitamins, minerals, and Detoxify's proprietary blend of cleansing herbs. Xxtra Clean is easy to use, in Tropical Fruit Flavor.

Detoxify Xxtra Clean Directions:

- Pick your day for using your Xxtra Clean detox drink for intensive cleansing.
- Begin your cleansing program with Xxtra Clean.
- Shake the Xxtra Clean well and drink entire contents of the bottle.
- Xxtra Clean begins to immediately work to help support your body's natural cleansing process.
- Wait 15 minutes. Refill the Xxtra Clean bottle with water — shake and drink.
- Urinate frequently. Frequent urination indicates that you are experiencing optimal cleansing with Xxtra Clean.
- Continue to drink plenty of water throughout the day to extend your cleansing program for hours.

To further enhance Xxtra Clean's effectiveness:

- Extend your intensive cleansing program with Xxtra Clean by taking PreCleanse Herbal capsules for the 2 days prior to using Xxtra Clean.
- Take Detoxify Constant Cleanse herbal capsules every day throughout the year. Periodic intensive cleansing is most effective when coupled with an ongoing detoxification routine.
- Eat light meals including fruits, vegetables, and fiber during your cleansing program with Xxtra Clean.
- Exercise regularly, and drink at least 12 8oz glasses of water in the days prior to using Xxtra Clean.

DETOXIFY MEGA CLEAN

http://www.detoxify.com/wp-content/uploads/2013/10/MegaClean_Box_39545_r05_80C_Bottle-1000sq.jpg

Detoxify Mega Clean Directions:

- Detoxify Mega Clean is our most powerful cleansing herbal detox Drink. Mega Clean is formulated for people with higher toxin levels, and limited time for their intensive cleansing program. The 32oz ready-to-drink formula is packed with vitamins, minerals, and Detoxify's proprietary blend of cleansing herbs. Mega Clean is easy to use, in Tropical Fruit Flavor
- Pick your day for using your Mega Clean detox drink for intensive cleansing.
- Begin your cleansing program with Mega Clean.
- Shake the Mega Clean well and drink entire contents of the bottle.
- Mega Clean begins to immediately work to help support your body's natural cleansing process.
- Wait 15 minutes. Refill the Mega Clean bottle with water — shake and drink.
- Urinate frequently. Frequent urination indicates that you are experiencing optimal cleansing with Mega Clean.
- Continue to drink plenty of water throughout the day to extend your cleansing program for hours.

To further enhance Mega Clean's effectiveness:

- Extend your intensive cleansing program with Mega Clean by taking PreCleanse Herbal capsules for the 2 days prior to using Mega Clean.
- Take Detoxify Constant Cleanse herbal capsules every day throughout the year. Periodic intensive cleansing is most effective when coupled with an ongoing detoxification routine.
- Eat light meals including fruits, vegetables, and fiber during your cleansing program with Mega Clean.
- Exercise regularly, and drink at least 12 8oz glasses of water in the days prior to using Mega Clean.

STINGER 7 day Permanent Cleanser

- Avoid toxins for 48 hours before taking the Stinger.
- Take one tablespoon in the morning and another tablespoon in the evening.
- If you have a "high toxin level" or are near or over 230 lbs double up on the 7 day making that 2 tablespoons in the morning and again at night.
- Make sure you are measuring with a kitchen Tablespoon, not just any spoon and that the measurements are accurate. You will have a couple tablespoons left over on the last dosage which you want to take all together on the last dose.
- This is a "permanent" detox. You will be clean and stay clean unless you put toxins into your body again.
- **This product will not work for those taking prescription medications, if possible, be off prescriptions for 48 hours before use just like you would toxins. "Prescriptions" also include anything over the counter such as Tylenol, Motrin, or cold medicine.**

Stinger Mouthwash:

- For best results use right before u go in, like in the parking lot, swish and swallow the entire bottle coat your mouth as best you can.
- Do not eat or drink anything after you use this product, it has a 3 hour window, but best time to use it is as close to the appointment time as possible...

Stinger Folli-Kleen Hair Shampoo:

- Wash your hair like you normally would.
- While your hair is damp apply Folli-Kleen and massage into your hair VERY WELL for 10 minutes straight.
- Let the Folli-Kleen soak into the hair for an additional 10 minutes for a total of 20 minutes.
- Rinse your hair after twenty minutes and you are good to go. Your hair is permanently clean.
- You can use any other products such as hair spray conditioner or cream after you are finished rinsing out the Folli-Kleen.

WHIZZINATOR

Preparing the Whizzinator for Wet Play:

- Check clip valve before filling to make sure it's closed
- Follow the included directions to mix the dehydrated urine
- In a cup, add about 80cc of tap water and mix the dehydrated urine.
- Fill your syringe with the appropriate amount of synthetic urine
- Unscrew the cap on the vinyl bag and fill the bag with your syringe
- Remove any air from the bag
- Screw back on the cap and make sure it is secure
- Open up the sticky heat pack
- Attach the heat pack to the vinyl bag on the opposite side of temperature strip
- Make sure it's attached to the bag, and not the belt.

Wearing your Whizzinator:

- Place the temperature strip side of the vinyl bag to your skin to help maintain heating
- If you did this correctly, the empty side of the vinyl bag should be against your skin and the heat pack should be on the opposite side of vinyl bag.
- The heating pad will keep the synthetic urine at a body temperature between 98 and 100 degrees for up to 8 hours
- Using the elastic straps you can adjust the fit of the Whizzinator to feel comfortable on you

Urinating with the Whizzinator:

- It's recommended that you practice using your Whizzinator with water before you actually use it
- Using one hand, hold the Whizzinator prosthesis and pinch it just below head to prevent any leaking or flow of urine, now using your other hand open the clip valve
- Once the clip valve is open, gravity will take over and allow the urine to flow naturally
- Do not use soaps or cleaning chemical to clean the bag, this will add new chemicals to any future Synthetic Urine you use and ruin the experience when you use it again

QUICK FIX

- Prior to leaving the house, open the top of your Quick Fix bottle and microwave for 10 seconds. If you can't microwave your Quick Fix synthetic urine, simply tape or rubber band the heating pad to the bottle at least one hour before you test.
- After heating, close the top and shake the Quick Fix bottle. The temperature strip should read between 94°-100°F. If the synthetic urine is below 94°F place in microwave for a few more seconds. If the temperature is above 100°F allow the bottle to cool.
- Tape or rubber band the heating pad to the Quick Fix bottle to maintain temperature until your drug test. When placed inside your clothing pocket, the heating pad will keep the Quick Fix synthetic urine near 100°F for at least 6 hours.
- Check that your Quick Fix synthetic urine is between 94°-100° F, shake and pour the contents of the Quick Fix container into the urine collection cup at the facility.
- QuickFix synthetic urine can be reheated with hand warmers an unlimited amount of times.
- Some settling may occur.
- Shake bottle before and after heating and right before you test the product.