

# Dometic Chill Box & LogTag Datalogger Kit

## For off-site vaccinations

The kit is made up and includes the following:

- 1x LogTag TRED30-16R temperature datalogger
- 1x LogTag ST-100K 1.5m probe with glycol bottle
- 1x LogTag mounting bracket
- 1x LogTag USB-HID Interface Cradle (if not an existing LogTag user)
- 1x Dometic Chill Box - 13L, 22L, 33L or 42L model



The kit comes completely installed and ready to use. Please read the instructions before use. The LogTag TRED30-16R datalogger is mounted on the outside of the Dometic chill box using a specially designed mounting bracket. The temperature probe is connected to the logger and is inserted through the top of the lid. The probe tip is inserted in a small plastic bottle which hold ethylene glycol. This reduces the response time of the probe during frequent lid openings and better represents the thermal kinetics of the product being stored. The LogTag data logger is configured and data downloaded by connecting to the USB-HID interface cradle and using specific LogTag Analyser software, which is available free from [www.logtagrecorders.com](http://www.logtagrecorders.com)

**The Temperature Shop**

P.O. Box 56-156 Dominion Road, Mt Eden, Auckland 1024, New Zealand  
Tel +64 9 265 0306 Email: [info@temperature.co.nz](mailto:info@temperature.co.nz) Web: [www.temperature.co.nz](http://www.temperature.co.nz)

 **BUTLER**  
TECHSENSE