


Anti-Rabbit Ig:DyLight® 488

Goat Polyclonal

Cat. # RS3261 **Size** 0.5 ml

Background

DyLight® fluorophores have absorption maxima covering the entire visible light spectrum (350 nm to 777 nm), as well several key near-infrared and infrared wavelengths. Both the absorption and emission properties of the DyLight® fluorophores match the excitation and detection wavelengths of common fluorescence instrumentation used for immunocytochemistry, immunohistochemistry, and flow cytometry. The DyLight® fluorophores exhibit higher fluorescence intensity and photostability than Alexa Fluor, CyDye and LI-COR Dyes in many applications and remain highly fluorescent over a broad pH range (pH 4-9). Additionally, the water solubility of the DyLight® fluorophores allows a high dye-to-protein ratio to be achieved without causing precipitation of conjugates. These favorable properties and high anisotropy value, as well as a high cross-section for two-photon excitation, make these fluorophores attractive as fluorescent probes in a variety of fluorescence methods.



γ-Catenin Immunocytochemical labeling of phosphorylation in A431 cells treated with pervanadate. The fixed cells were labeled with rabbit polyclonal anti-v-Catenin (Tyr-550) phosho-specific and the antibody was detected using Goat anti-Rabbit secondary antibody conjugated to DyLight® 488.

Background References

Sarkar P. et al. (2010) J Photochem Photobiol B. 98(1):35. Alkhatatbeh MJ (2011) J Thromb Haemost. 9(4):844.

Applications Species Reactivity Specificity

IHC 1:200

End user should determine optimal dilution for their particular applications and

Western blot membranes were incubated with diluted antibody in 5% non-fat milk, PBS, 0.04% Tween20 for 1 hour at room temperature.

DvLight® is a trademark of Thermo Scientific. Inc. and its subsidiaries.

This antibody has been pre-adsorbed with various immunoglobulins from nonrabbit species before affinity-purification using rabbit IgG coupled to agarose beads. Purified goat polyclonal antibody was conjugated to DyLight® 488 fluorophore (Excitation = 493 nm; Emission = 518 nm). This secondary reagent can be used to detect rabbit immunoglobulins in various immunofluorescence applications, such as immunocytochemistry, immunohistochemistry, and flow cytometry.

*All molecular weights (MW) are confirmed by comparison to Bio-Rad Rainbow Markers and to western blot mobilities of known proteins with similar MW.

Immunogen

The DyLight® 488-conjugated goat polyclonal secondary reagent reacts with the heavy chains of rabbit IgG, and the light chains of most rabbit immunoglobulins. The secondary reagent has minimal cross-reactivity with goat, horse, human, mouse, rat, and swine immunoglobulins.

Buffer and Storage

Goat polyclonal antibody is supplied in phosphate-buffered saline, 50% glycerol, 1 mg/ml BSA and 0.05% sodium azide. Store at -20° C. Stable for 1 year.

Related Products

PF7501 Phalloidin:FITC Reagent PF7551 Phalloidin:TRITC Reagent

MS3011 Anti-Mouse Ig:DyLight® 488 Goat Polyclonal MS3031 Anti-Mouse Ig:DyLight® 594 Goat Polyclonal RS3271 Anti-Rabbit Ig:DyLight® 594 Goat Polyclonal

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE

www.ecmbio.com toll-free: 1-800-859-8202 telephone: 859-879-2075 info@ecmbio.com