

## Product Description

A pull door handle that automatically applies sanitising hand gel onto the user's hand when they open the door

## Operation

When a person takes hold of the Purehold PRO to open the door, sanitising hand gel is automatically released onto their hand through the 19 individual release valves positioned on the handle surface

## Model

PH-HH-001 – Purehold PRO

## Specifications

- Intelligent gel delivery system (19 release valves)
- Secure anti-theft gel storage
- Robust aluminium handle with strong stainless steel frame (compliant with BS 8424)
- Large gel capacity - each refill lasts around 2,000 "pulls" of the handle
- Low gel level indicator (red flotation ball)
- Easy replacement of disposable hand gel bottles
- Low maintenance design
- Approved for use on FD30 and FD60 timber and steel fire doors (tested to BSEN1634-1: 2008\* and BSEN1363-1:1999)

## Benefits

- Proven to significantly increase hand hygiene compliance and clean hands by 87%
- 'Intelligent' delivery system ensures that the volume of alcohol gel delivered is always appropriate to the size of the user's hand
- Proven to be 98.5% cleaner than standard door handles
- Uses virucidal gel for superior performance
- Can be retrofitted onto doors using original pull handle fixing points (mounting hardware included)

## Installation

Undertaken by a trained Pure Hold Installation Engineer  
In the majority of installations, the Purehold PRO can be retrofitted onto doors without drilling holes or screwing into the doors

## Warranty

12 months from date of installation

## Gel Replacement

Extremely fast and simple to replace the Pure Hold Sanitising Hand Gel Disposable Bottle. When the red floatation ball becomes fully visible in the reservoir window it indicates that the gel bottle is empty and that it needs to be replaced

## Pure Hold Sanitising Hand Gel

The Pure Hold Sanitising Hand Gel comes in a 850ml disposable bottle.  
PH-GL-001 — 850ml Pure Hold Sanitising Hand Gel  
PH-GL-009 — Pack of 9 x 850ml Pure Hold Sanitising Hand Gel

## Order Information

To order simply: call us on +44 (0)2392 299100 e-mail us [info@purehold.co.uk](mailto:info@purehold.co.uk) or visit our website [www.purehold.co.uk](http://www.purehold.co.uk)




\* satisfies the provisions of BS 476: Part 22: 1987 until harmonised product standards are published

## Dimensions

Height - 559mm

Width - 91mm

Depth - 97mm



## Table of Components / Replacement Parts

Sub-Part Numbers	Sub-Part Name	Description
PH-RP-KEY	Security Key	Security key that allows user to open the 1) Top Lid (b) and 2) Bottom Cover (k)
PH-RP-055	Drip Tray	Collects any excess gel that may run down the Handle during high usage. Remove this part to gain access to 1) Excess gel that may have been captured and 2) the Non-Gel Handle fixing screw
PH-GL-009	Gel Refill Bottles	Disposable gel refill bottle (sold separately)

In order to help keep the long term support costs for the handles low, individual replacements parts are available. Items shown above with sub-part numbers can be purchased as spares.