

PROTOCOLO DE APLICACIÓN

TÉCNICAS DE DERMATOLOGÍA ESTÉTICA PARA UN RESULTADO REDUCTOR PROFESIONAL

Protocolo de aplicación desarrollado en colaboración con el Centro Médico-estético Kursal (www.kursal.es), para optimizar y aumentar los efectos reductores y anticelulíticos del innovador sistema liporeductor Glo910.

El protocolo de aplicación permite reproducir fácilmente maniobras análogas a las realizadas en Centros especializados y conseguir resultados profesionales en tu propia casa.

PROGRAMA DE EXFOLIACIÓN CORPORAL: PIEL RENOVADA Y LUMINOSA

El cepillado en seco tiene muchos beneficios para la piel, el órgano más extenso del cuerpo, este tratamiento exfoliante activa la circulación y el metabolismo, a la vez que deja la piel libre de células muertas, más suave y tersa.

Con glo910 la exfoliación se realiza en **menos tiempo y sin esfuerzo** que un tratamiento manual, arrastrando las células muertas y las impurezas, a la vez que tonifica el cuerpo, activa el sistema linfático y la microcirculación sanguínea, relaja, detoxifica y aumenta la eficacia del tratamiento anticelulítico.



REALIZAR EL TRATAMIENTO:

- Con el cabezal **BODY BRUSH**
- Utilizar el cabezal **LIFTING-CLEAN** para exfoliar la zona del cuello y escote
- Sin cremas, sobre la piel limpia y seca.
- Intensidad del masaje: Regular la intensidad de la exfoliación en función del estado de la piel, zona del cuerpo y sensibilidad particular.
- Tecnología **SLIM Tech®**: Opcional

Siempre hay que cepillar la piel en dirección al corazón, empezando por los pies y subiendo hacia la cadera en el caso de las piernas, en brazos de la muñeca al hombro. Insistir en las zonas más rugosas y ser cuidadoso con las zonas delicadas. Regular la intensidad de la exfoliación en función del estado de la piel, zona del cuerpo y sensibilidad particular.




Evitar dañar la piel por una excesiva exfoliación (1 o 2 veces por semana son suficientes). Realizaremos un masaje muy suave durante 5 minutos por toda la zona.



Después, es aconsejable ducharse y aplicarse agua fría para finalizar la ducha, al menos en las piernas. Terminar aplicando una crema o aceite hidratante.



Además, este tratamiento deja la piel con la base perfecta para aplicar cremas hidratantes o bronceadoras, ayudando a estos productos a lograr el acabado perfecto.



APPLICATION PROTOCOL

AESTHETIC DERMATOLOGY TECHNIQUES FOR PROFESSIONAL FAT REDUCTION RESULTS

Application protocol developed in cooperation with Kursal Medical-Beauty Center (www.kursal.es) to enhance fat reduction and anticellulite effects of the innovative glo10 device.

This application protocol helps reproducing similar way of massage procedures done at Specialized Centers achieving professional results on the comfort of your home.

BODY EXFOLIATION PROGRAM

Dry brushing has many benefits for the skin, the body's largest organ. This exfoliating treatment stimulates circulation and the metabolism whilst leaving the skin free from dead cells, softer and smoother.

Glo910 achieves professional exfoliation results, faster and with less effort than manual treatments. Exfoliation removes dead cells, tones the body, activates lymphatic system and blood circulation, while relaxing and detoxifying, enhancing anticellulite treatment efficacy.



PERFORM THE TREATMENT:

- Using the **BODY BRUSH** head
- Use **LIFTING-CLEAN** for neck and neckline exfoliation
- Without creams, on dry clean skin
- Massage intensity: Adjust the intensity of exfoliation based on the condition of the skin, the body area and particular sensitivity.
- **SLIM Tech®** technology: Optional

The skin should always be brushed toward the heart, starting with the feet and moving up to the hips in the case of legs, and from wrist to shoulder for arms. Persist on rougher areas and be careful with sensitive areas. Adjust the intensity of exfoliation based on the condition of the skin, the body area and particular sensitivity.



Avoid damaging the skin through excessive exfoliation (once or twice per week is sufficient). Gently massage the entire area for 5 minutes.



Afterwards, it is advisable to shower and finish the shower with cold water, at least on the legs. Finish by applying a moisturizing cream or oil.



This treatment brings your skin to perfect conditions to use moisturizer and suntan creams, enhancing their results.

