



MOBILE OPERATING CHAIR



The OS1 is designed for mobility and flexibility in outpatient surgery.

This new chair meets the patient positioning and mobilization challenges posed by the growing number of outpatient surgery centers. It is ideal for surgical preparation, transporting patients to and from the operating room, positioning after surgery, and deployment in the recovery room after anesthesia.

Ergonomics, functionality, mobility and patient comfort form a unique combination of chair and operating bed. Complete ease of access and simple disinfection ensure maximum hygiene.

The **OS1** can also be installed in small spaces where conventional beds will not fit.

JEDMED
Employee-Owned

FEATURES

Product Safety:

OS1 is independently tested and certified

Ergonomics:

Fixed and retractable armrests for simple repositioning before and after surgery

Minimal Space Requirement:

In a seated position with the footrest retracted, the OS1 has a minimal footprint

Mobility:

5" casters provide comfortable and effortless mobility

Choice of Headrests:

- Ergonomic Full Headrest
- Full Spherical
- Half Spherical



Flexibility:

Infusion rod holder with right/left swiveling position - one handed operation

Hygiene:

Compliant with stringent cleaning regulations

Seated and Reclining Comfort:

Comfortable, 24" wide seating for pre- and post-operative waiting and recovery

Stability:

High weight capacity up to 500lbs. Safe handling during patient transport and for lateral repositioning

Fixation:

One pedal for central braking system

COMFORT, FUNCTIONALITY AND ERGONOMICS

With stepless motorized adjustment and intelligent, intuitive operation via the control pad, the OS1 Mobile Operating Chair delivers outstanding seated and reclining comfort in all situations.

Seated Position



Relaxed Position



Flat/Bed Position



Trendelenburg Position



COMFORT FOR OUTPATIENT SURGERY


The OS1 gives you the opportunity to individualize your chair according to your precise needs.



Ergo headrest with vario double joint



Full spherical cap for quick assist




Half spherical cap for quick assist



Change headrest pad with quick assist coupling



Any position can be fixed with a vario double joint



Armrests fully retractable



Standard bars for clinical accessories

