

# Stand-to-Pee (STP) Packers and Devices



## Available at Trans Tool Shed

2427 N Murray Avenue, Milwaukee, Wisconsin  
[transtoolshed.com](http://transtoolshed.com)



### Pstyle \$12

- Hard plastic
- Easy to use while fully clothed
- Difficult to conceal in pocket
- Comes in seven colors



### Go Girl \$10

- Bendable, foldable silicone
- Difficult to use while fully clothed
- Easy to fold and conceal in a pocket.
- Comes in pink or camo (tan)



### The Fenis \$15

- Bendable, foldable silicone
- Somewhat difficult to use while fully clothed
- Easy to fold and conceal in a pocket.
- Comes in numerous colors and 4 flesh tones (vanilla, caramel, olive, chocolate)



### Mr. Fenis \$25

- Bendable, foldable silicone
- Larger than Fenis and with added head
- Somewhat difficult to use while fully clothed
- Easy to fold and conceal in a pocket.
- Comes in 4 flesh tones (vanilla, caramel, olive, chocolate)



### Number One \$52

- Flexible silicone
- Funnel is built into the design and a hole runs through the length of the shaft; no tubing
- Doesn't wear/pack like a typical packer due to shape
- 3 colors: vanilla, caramel, and chocolate



### **STP Fitz 2.0 Mr. Sofite \$38 & \$39**

- Mr. Softie soft packer retrofitted with a tubing/funnel combo through the shaft
- Tubing and funnel are made of silicone
- Tubing/funnel is a single piece, eliminating glue and preventing breakage
- Funnel is soft enough to collapse into itself partially during regular packing
- Available in small and large
- 2 colors: vanilla and caramel



### **Packer Gear STP \$15**

- Somewhat firm 100% silicone, but entirely hollow so it can be folded into packing position
- Funnel is built into the design; no tubing
- Not as "squishy" as most soft packers, but firmness of funnel may help some with use
- Hollow core may help those who experience overflow with other designs
- 2 colors: vanilla and caramel



### **Packer Gear Ultra Soft Silicone STP \$35**

- Flexible silicone
- Funnel is built into the design and a hole runs through the length of the shaft; no tubing
- Entirely hollow; may work for some with metaoidioplasty, depending on size/length
- 4 colors: ivory, tan, brown, and dark brown



### **Sam the STP \$75**

- Flexible silicone
- Funnel is built into the design and a hole runs through the length of the shaft; no tubing
- Optional silicone insert can be purchased for penetration function
- 4 colors: cashew, caramel, hazelnut, and chocolate



## Model A \$80

- Flexible silicone
- Available in circumcised and uncircumcised
- Largest of the "Model D" family of products
- Funnel built into the design, a hole runs through the shaft; there is no tubing
- 3D testicles
- 3 colors: vanilla, caramel, and chocolate



## Model D \$75

- Flexible silicone
- Available in circumcised and uncircumcised
- Mid-size of the "Model D" family of products
- Funnel built into the design, a hole runs through the shaft; there is no tubing
- 3D testicles
- 3 colors: vanilla, caramel, and chocolate



## Sport \$70

- Flexible silicone
- Available in circumcised and uncircumcised
- Smallest of the "Model D" family of products
- Funnel built into the design, a hole runs through the shaft; there is no tubing
- 3D testicles
- 3 colors: vanilla, caramel, and chocolate



## Model E (youth size) \$50

- Flexible silicone
- Available in circumcised and uncircumcised
- Designed specifically for trans youth, for whom an adult-size STP is too large
- Funnel built into the design, a hole runs through the shaft; there is no tubing
- 3D testicles
- 3 colors: vanilla, caramel, and chocolate

# **What to look for when choosing an STP device**

Much of the advice below is geared toward trans men, who often have specific needs when it comes to privacy and safety when using an STP. But some advice will still be relevant to all users.

## **1. Will I be using the STP in privacy, or in a more public place like a urinal, campsite, or worksite?**

If you typically plan to use your STP in a bathroom stall, or at home, or any other mostly private place, you might not care so much about how your STP looks, but rather how well it functions. However, if you are planning on using your STP in a men's room at a urinal, or in a semi-public setting such as a campsite, a boat, or a worksite, you may be concerned about what others will see if they catch a glimpse of your STP. If that is the case, you might want to consider an STP packer, or at the very least something that appears to be flesh-toned and that is easily tucked into your pants after use.

## **2. Do I want an STP that looks like a penis?**

As noted above, certain situations (like urinal use or other semi-public use of your STP) might be safer and less stressful with an STP that looks like a penis. But even for those who always pee in privacy, having an STP that looks like a penis may help them feel more at home in their bodies. Many trans guys choose prosthetic penises with an STP option for this added sense of well-being.

## **3. How will I clean and store the STP when not in use?**

Some STPs are designed to be easily folded and hidden in a pocket or a bag, while others can be a bit more unwieldy. Keep in mind how you plan to get the STP to the bathroom with you, and whether you'd like to rinse it out and store it again immediately after use.

## **TAKE NOTE: MOST STPs TAKE PRACTICE TO MASTER!**

All STPs, especially STP packers, will take some practice to get the hang of it. It is best to first practice at home, without your clothes on in the shower to get used to correct placement and flow. While still at home, practice in front of the toilet without clothes, and then finally try it at in front of the toilet while clothed. This is crucial if you plan on using a public urinal! It is better to have any accidents or mishaps at home than when you are out in public.