

ALPHASONIC II

MASSAGE CHAIR

< timeto**relax**



 braintronics®
brainwave stimulation system


CASADA
HEALTHCARE

SPECIALS



BRAINTRONICS®

Braintronics® technology synchronises your brainwaves through audio stimulation. Combined with a specially developed massage programme your stress-ridden body will be able to relax.



SHIATSU

Shiatsu (finger press) is a form of body therapy originating in Japan. The massage technique consists of soft, rhythmic, far-reaching stretches and rotations. The aim of Shiatsu is to simulate individual parts of the body and to mobilise the muscles.



2D MASSAGE

The intelligent 2D massage mechanism allows the massage chair to measure your back to determine a more exact position of the sections of the spine.



KNEADING

Skin and muscles are either clasped and kneaded between the thumb and the index finger or with both hands. This massage technique is used in particular to ease tension.



FLAP

Circulation in the muscles is stimulated and the tendons are encouraged to stretch by means of light tapping with a cupped hand, which also reduces tension at the same time.



TAPPING

Short, tapping movements are performed with the edge of the hand, the palm or the fist. This promotes blood circulation in the skin and softens tense muscles. If tapping is performed at lung height, this can improve the release of mucus in the lung.



KNEAD & FLAP

A massage combination of kneading and flapping.



ROLLING

Pleasant alternation between relaxation and releasing the muscles – also perfectly suited for a final massage, for relaxation and recovery of the muscles.



CARBON HEATING

Deeply penetrating Carbon infrared heat, blood vessels expand, the blood flow increases and the muscle tone is reduced. In addition, it carries out a highly beneficial effect on nerve cells, which are responsible for reporting pain to the brain.



GUA SHA

The scraping technique of -GuaSha- with gentle pressure opens up superficial blood vessels and stimulates blood circulation. GuaSha is often applied in combination with acupressure 'shu points' for additional stimulation of the whole body for an extra soothing massage experience.



AIR COMPRESSION

Rising and receding airbags have a pumping motion which produces natural muscle reflexes. This is particularly beneficial in the calf and foot area.



L-SHAPE

You can expect a profound full body massage from the neck down to the buttocks.



S-LINE

With this automatic program the massage heads are working very close to your body, they massage the whole back down to your buttocks. The S-shaped mechanic is constructed after the natural anatomy of the spine.



ZERO SPACE

For the lying-position the massage chair moves forward and therefore requires very little distance to the wall behind. Space saving and comfortable!



ZERO GRAVITY

You put yourself in the position of an astronaut – that "weightless feeling".



ADJUSTABLE PILLOW

If necessary, the double cushion can be placed to reduce the neck and shoulder massage intensity.



AUTOMATIC

9 specifically developed automatic programmes for the whole body, which offer a big range of different massage techniques. Additionally you can also turn on the heating function, if you'd like.



MEMORY

Using the intelligent memory function the massage chair remembers your favourite massage exactly and begins immediately without a scan. You have the opportunity to save two programmes with different buttons (M1 & M2).



BACK AREA

An individual applicable to any size massage area.



ADJUSTABLE FOOTREST

The footrest can be extended, this way it is adaptable to different body sizes.



RECLINING POSITION

The adjustability of the backrest and footrest contributes to personal well-being of the body position.



ADJUST SHOULDER

The massage heads can be moved upward or downward to a suitable position in the shoulder area.



AIR INTENSITY

Adjustable airbag massage intensity.



AIR AREA

Three air pressure modes (full body, arm and leg / lower body).



WAIST STRETCH

Airbags in the lower area of the back ensure gentle stretching of the pelvis. Effectively counteracts tensions.



SPOT

The massage heads can massage a particular spot.



SPEED

Massage speed is adjustable.



WIDTH

The distance between the both massage heads can be set (narrow, medium and wide).



ACUPRESSURE POINTS

Key acupressure points are massaged.



MUSIC

Further relaxation and fun with the music function.



BLUETOOTH

Wireless data transmission.



HUMAN HANDS FEELING

Specially designed massage heads imitate the hand motion of a professional masseur. The maximum stretch space is 6,5 cm, and the maximum stretch angle is 41 degrees, so the mechanical hands can do effective massage to human cervical vertebra thoracic vertebra and lumbar vertebra.



O2 IONIZER

Polluted air around you will be neutralized and you will breathe purified oxygen during your massage.



USB-CHARGER

While you relax, the integrated USB port can supply, for example, your phone with power.



QUICK START

With one touch the quick start buttons in your armrest allow you to control the most important functions.



YOGA

Targeted stretch massages activate the body and have an invigorating effect on the whole of the body.



REFLEXOLOGY

Reflexology regenerates foot muscles and has a harmonizing effect on almost all organs inside your body.



braintronics®

brainwave stimulation system

The human brain is a complex organ consisting of millions of little nerve cells, called neurons. Each cell constantly produces tiny electrical impulses. Every thought, every sound and every impression, everything is transmitted and communicated by those small electrical impulses. Electricity is, so to say, the language of the brain. The sum of those impulses can be detected and diagrammed, e.g. by using an EEG. Charted on a graph this electrical activity forms a wave-like pattern, known as brainwaves.

With braintronics®, combining audio stimulation and guided meditations with a soothing massage, these brainwaves can be stimulated in order to reach different states of relaxation in the brain. Whether you want to sustainably enhance your concentration, or just relax and pamper yourself – everything is possible with a massage chair from the Casada family.



BODY + MIND



braintronics[®]
brainwave stimulation system

ALPHASONIC II



iPhone



Android



ALPHASONIC II




The next step towards perfection – **AlphaSonic II from Casada**

The AlphaSonic II combines state-of-the-art technology with tried-and-tested/proven quality and offers a completely new way to experience relaxation. The new premium massage chair from the Casada family impresses with a larger seat and a striking, creative design. Durable and environmentally friendly microfiber PU leather does not only feel good, but also facilitates cleaning and ensures exceptional seating

comfort. With the quick start panel located in the armrest, all basic functions are easily controllable, and you can simply put the remote control aside. The new braintronics® App makes holistic relaxation easier and more accessible than ever. Dream away, take a break from the life's daily hassles and allow body and mind to relax.



CLASSY DESIGN AND WONDERFUL RELAXATION



Headphones
additionally
available



BRAINTRONICS



TECHNICAL DATA

Dimensions upright	81 x 119 x 147 cm
Dimensions reclining	81 x 92 x 176 cm
Weight	96 kg
Voltage	220-240 V ~ 50/60 Hz
Nominal power	150 W
Rated runtime	20 minutes
braintronics®	21 minutes
Certificates	

AVAILABLE COLOURS



Width of seating area 50 cm | Width of shoulder area (inside) 60 cm



Casada International GmbH
Obermeiers Feld 3
33104 Paderborn | Germany
www.casada.com