

## **PRODUCT OVERVIEW**



TGA Listed Hospital Grade Surface Disinfectant that creates a protective barrier with 28-day residual activity. Extensively tested to be effective against COVID-19 and 99.99% of germs. Can be used as a hospital grade disinfectant, or diluted for regular, daily use.

### **PRODUCT BENEFITS**



The long-lasting barrier reduces cleaning time and cleaning costs and has been tested to be effective against 99.99% of germs including COVID-19. Remains on surfaces for more than 60 scrub cycles with a sponge and 40 scrub cycles with microfibre cloths.

## **WHERE TO USE**



Any surfaces that require cleaning and sanitising.

# **TECHNICAL SPECIFICATIONS**

COLOUR: CLEAR pH: 3.0>5.0 ODOUR: NO SCENT

## **HOW TO DILUTE**



For TGA Listed (Hospital Grade) Disinfectant: Do Not Dilute

**For Standard Disinfectant Cleaner:** Add 150ml to a 750ml spray bottle and top up with clean water.

## **DIRECTIONS FOR USE**



#### For TGA Listed (Hospital Grade) Disinfectant:

- Spray, mist or wipe directly onto surface and allow to dry for long lasting protection.
- For 28-day residual activity, reapply after first application has dried.



#### For Standard Disinfectant Cleaner:

• Use as a daily cleaner/disinfectant to maintain the antimicrobial coating.

## **PPE & SAFETY**



**PPE:** Wear protective gloves & eye protection.

**Safety:** Always take care when handling any chemicals.

Refer to the SDS for additional information at www.fogeliminator.co.nz



