TANITA® BC 541

It's How We Measure.

9-in-1 Body Composition Monitor

Equipped with 9 basic body composition features, The BC-541 allows you to have a quick overview of your body composition situation at a glance.

The BC-541 is easy to use, users can check each body composition reading simply by clicking the corresponding feature buttons on the monitor. It also features a durable dark grey platform robust for frequent use.

Convenient Data Collection

| \bigcirc | |
|------------|--|
| | |

Your body composition measurements are automatically saved between the 4 available user profiles.



Track your progress with easy-to-read results for helpful assessments.

Set your own target and achieve your personal best.

TANTA TA

RRP \$199

| Weight | Image: A second s |
|--------------------------------|--|
| Body Fat % | Image: A second s |
| Total Body Water % | Image: A second s |
| Muscle Mass | Image: A second s |
| Physique Rating | \sim |
| Bone Mineral Mass | Image: A second s |
| Basal Metabolic Rate (BMR) | |
| Metabolic Age (Up to 50 years) | |
| Visceral Fat | \sim |
| | |

Features

| Recall Function |
|---|
| Guest Mode |
| Athlete Mode |
| 4 User Memory |
| Children's Body Fat Range (7-17 years) |

>>>>>

Product Specifications

| Weighing Capacity: | 150kg |
|--------------------|--------------|
| Weight Increments: | 0.1kg |
| Dimensions: | 302x307x43mm |
| Product Weight: | 2.4kg |
| Warranty: | 5 years |

Want to keep track of results for even more users via Guest Mode?

Download the Healthy Range Charts to keep track of your health.

https://tanitaaustralia.com/pages/download-healthy-range-charts

Can't find what you're looking for?

Visit the Tanita Australia website: www.tanitaaustralia.com

