

I N O G E N O N E^{. G5}

The Inogen One G5, with flow settings from 1-6, is designed for most supplemental oxygen users 24/7. Advanced diagnostics take the guesswork out of patient support, as the Inogen One G5 is bluetooth enabled and compatible with Inogen Connect, the free mobile app. It can be charged at home or on the go.

Weight	4.7 lbs. (includes single battery)
Size	Length: 7.19" Width: 3.26" Height: 8.15" (includes single battery)
Oxygen Flow	Pulse dose delivery system Six flow settings 1, 2, 3,4,5,6
Power	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): Up to 6.5 hours [△] Recharge up to 3 hours with AC or DC power Duration (double battery): Up to 13 hours [△] Recharge up to 6 hours with AC or DC power
Noise Level	38 dBA*
Operation	Simple control functions and easy-to-read LCD display Inogen Connect App
Use	Designed for 24/7 use and meets FAA acceptance criteria for use on board aircraft.**

* At flow setting 2

** Based on technical performance test data.

△ At Flow Setting 1

Manufactured by:





PM230196 v12303