



## Plating Limit For Cyanideless 18KT Rose Gold Plating Solution

CLRGPS is a ready-to-use plating bath at the concentration of 2 g/l of gold. No preparation is required.

The initial gold concentration; for example, with a bath operating at a concentration of 2 g / l it can be used upto a maximum consumption of **0.4 g / l** of Gold.

<b>Deposit Data</b>	
Purity (%)	75
Hardness [HV 0.01]	160 - 220
Density [g/cm <sup>3</sup> ]	19
Thickness from-to [μm]	0.02 – 0.20
Aspect	Shiny
Color	18 KT Pink