



CANINE IgG SOLUTION

TDS-CA60

Effective Date: February 27, 2019

Rev. 000

Page 1 of 1

METHOD OF PURIFICATION: Equitech-Bio's purified immunoglobulin (IgG) preparations are isolated from pooled normal serum by a proprietary fractionation method.

SPECIFICATIONS:

Form:	Liquid
Source:	Normal canine serum
Filtration:	0.2 μ absolute
Purity (Electrophoresis):	$\geq 97\%$
pH:	7.2 \pm 0.2
Buffer:	10mM sodium phosphate, 150mM sodium chloride pH 7.2
Concentration:	10mg/ml is the approximate normal concentration. Other concentrations up to 40mg/ml are available on request.
Preservative:	None

CATALOG NUMBER

CA60-0100

CA60-0500

CA60-1000

SIZE

100mg

500mg

1gm

SHIPPING CONDITIONS: <0°C

STABILITY AND STORAGE:

This sterile solution is stable for ≥ 3 years when stored unopened at $-20^{\circ}\text{C} \pm 5^{\circ}\text{C}$.

The material for this product was processed in the USA from the serum of healthy animals that had not been inoculated with or exposed to any infectious agents. During processing this material was subjected to a pH of ≤ 5 for at least 3 hours. At no time during collection or processing was the material commingled with any other material of animal origin.

This product is not to be used for direct or indirect exposure or inoculation into laboratory and domestic animals.

This product is sold for use in IVD, research, further manufacturing, or technical purposes only, and is not intended for direct use in diagnostic or therapeutic applications.

All of the dogs are vaccinated against rabies, canine parvovirus, canine distemper virus, and canine viral hepatitis virus.