
HOW TO USE 3 OUT

3 Hrs Out contains electrolytes and several innovative ingredients to help the animal metabolize energies to give extra muscle shape and definition.

GIVE THE ENTIRE 15CC TUBE OF 3 HRS OUT A FULL 3 HRS BEFORE THE ANIMAL WALKS INTO THE SHOW RING.

- Give the entire tube regardless of the size/weight of your animal.
- It takes an entire 3 hours to work.
- If you are guessing when you need to administer the product, guess early, NOT late.
- Once 3 Hrs Out begins to work, the animal will hold that look for 3 hours.
- If you show again outside of the 3 hour window, administer another tube so the animal will maintain that look.





HOW TO USE 3 OUT

FAQ'S

- Will using 3Hrs Out get me in trouble?

 3Hrs Out contains no drugs or banned substances and will not get you tested out at any show in the USA.
- How long does it take to see results of 3Hrs Out and how long will they last?

 It takes 3hrs for the product to reach full effect and the animal will hold this look for 3hrs.
- What should I change in my show prep when using 3Hrs Out?
 You do not need to change anything in your show prep
 - regimen in order to use 3Hrs Out. You do not need to increase fluid intake, add electrolytes, increase sugars, nothing.
- What if my show doesn't allow pastes?
 No problem! Find our recipe for Banner Bites at
 www.StockShowSecrets.com These need to be made at home,
 and you feed them to your animal at the show.

