

## **TRYPsin WASH (0.25%)**

Trypsin Wash (0.25%) is a 0.25% solution (Hank's Balance Salts) of Trypsin (1:250) without either  $\text{CaCl}_2$  or  $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$ . As specified by the embryo import protocols of many countries, trypsin is used prior to cryopreservation to remove/deactivate several different viruses that may be on the surface of the zona pellucida. This medium is packaged in 5ml biofreeze vials. Shelf life is 9 months from date of manufacture and the product must be stored frozen.

**Recommended use:** Embryos for export should be washed in Trypsin Wash (0.25%) in accordance with the regulations set by the importing country.

**Trypsin Wash (0.25%) should be stored frozen at all times.** Thaw immediately prior to use. **Do not refreeze.**



1615 NE Eastgate Blvd, Sec H  
Pullman, WA 99163 USA  
Phone 509.592-8144  
[www.abt360llc.com](http://www.abt360llc.com)