

## BIOFA Brush Cleaner

**solvent free**  
**Art. Nr. 0600**

### Properties

BIOFA Brush Cleaner is a solvent-free, bio-degradable cleaning solution for cleaning paint and lacquer instruments such as brushes, rollers, etc. For all natural resin paints, lacquers, glazers, oils and waxes. Cleans thoroughly but gently. Brushes stay soft and pliant.

### Contents

Water, potash saponified coconut and soy fatty acids, 5-15%, lavender oil (contains linalool).

### Use

Remove excess paint, lacquer, etc. immediately after use by wiping off on the edge of the container or striking it off. Place the tool in undiluted Brush Cleaner and let it soak for 5-10 minutes. Then wash out thoroughly and follow this by rinsing with clean water. Repeat as necessary. Heavily soiled items should be soaked in Brush Cleaner for 24 hours. BIOFA Brush Cleaner can be re-used. Leave the used cleaner in a container until the dissolved paint residue settles out. Then pour the clear solution on top into a clean container.

For very stubborn soiling that cannot be removed with Brush Cleaner, use BIOFA Thinner 0500.

### Storage

Store in a cool, dry and well-secured place. Protect from frost.

### Container

1 l PE- or PP-container

### Safety Instructions

Keep out of reach of children. Avoid eye contact. In case of eye contact, flush eye thoroughly and immediately with water. If ingested, seek medical advice immediately and show package or label to medical personnel.

### Disposal

Dispose of product residues and empty containers in compliance with local regulations.

Waste Incineration Ordinance [AVV]-  
Waste Code under European Waste Catalogue: 08 01 16

Die Angaben und Hinweise des Technischen Merkblattes sind verbindlich. Falls von diesen Vorschriften abgewichen werden muss, ist mit der anwendungstechnischen Abteilung der Fa. BIOFA vorher Rücksprache zu halten. Die allgemeinen Regeln der Bautechnik müssen eingehalten werden. Mit Erscheinen dieses Merkblattes verlieren alle bisherigen Angaben ihre Gültigkeit.