

Baseline AcuAngle Inclinometer

Product Information



Weight 0.132 kg

Dimensions 10.16 x 2.54 x 10.16 cm

Brand Baseline

Item No. W54668 [1013982]

On Sale

Description

Baseline AcuAngle inclinometer's pointer is damped by fluid to assure accurate range-of-motion measurements. Place inclinometer near joint to be measured; turn dial to zero; take joint through its range; read range of motion from dial. Adjustable legs (side-to-side) adapt to body contours.

https://www.3bscientific.com/ph/baseline-acuangle-inclinometer,p 906 17032.html