PLANTAR FASCIITIS PROGRAM

BRACE LAYER SYSTEM™

The

Comprehensive Solutions from the #1 selling brand for Plantar Fasciitis!

| | | 0 | |
|--|-----------------|---|--|
| Available in 5 Colors! | Usage Type | Active | The patented OS1 st FS6 Compression Foot Sleeves create an orthotic effect on feet, so when wearing an FS6 on 1 foot and not the other, you may experience a change in your gait and create a biomechanical imbalance. While this may not be true of other so-called foot sleeves or arch supports, our Compression Zone Technology® used to create OS1 st , was designed for this very specific reason. The evidence of this change was confirmed through foot specialists including Podiatrists, Certified Pedorthists & orthopedic surgeons, as well as millions of active people worldwide. |
| | Purpose | Performance & Prevention | |
| | Condition Level | Moderate and/or Recurring PF | |
| FESCIE Performance Foot Sleeve SINGLES | Usage Type | Recovery Resting / Light-Use | The Original FS6 Single Sleeve is now back as part of a complete Plantar Fasciitis treatment offering. The single is meant to be worn when you are in a resting position and/or recovering. The single FS6 can be a simple, less expensive treatment for mild cases or repeat onset of pain and swelling. It is the perfect treatment for the end of the day or after activity or excercise for relief and recovery. |
| | Purpose | Rehab & Recovery | |
| Available in Black Only | Condition Level | Mild - Moderate | |
| DS 55 Decompression Sleeve | Usage Type | Resting Sleep / Overnight | New overnight/resting PF treatment device unlike anything currently on the market! Recommended for more severe PF symptoms when chronic pain is present. This device will offer relief, especially in the morning after |
| COMING SOON! | Purpose | Pain Relief & Therapy | sleep. The DS6 Decompression Foot Sleeve is created using patent pending K-Zone Technology. This creates a comfortable, wearable overnight/resting therapy that offers benefits consistent with clinical taping. The DS6 offers an easier to use, more comfortable alternative to typical soft splints or socks. Therapists may choose to use the DS6 where taping was preferred and in some cases, instead of night splints. |
| Available in White Only | Condition Level | Moderate - Severe PF | |
| Available in 5 Colors! | Usage Type | Active | Combined benefits of the FS6 & CS6 Performance Sleeves relieves and prevents Plantar Fasciitis, eases heel pain, improves circulation, relieves shin splints, cramping & aids muscle recovery. Innovative design resulting in an easy-to-use orthopedic support with a thin profile to fit under any apparel. The perfect alternative to purchasing multiple pairs of compression socks so you can wear your favorite of socks over the sleeves. |
| | Purpose | Prevention Circulation | |
| | Condition Level | Mild - Moderate PF Before, During or After Activity | |