

HAINES® DISPOSABLE CURTAINS

Significant cost savings and improved infection control

Haines®
MEDICAL AUSTRALIA

NOW
AVAILABLE IN 6
STANDARD & 12
SPECIAL ORDER
COLOURS



Visual Appeal

Curtains always look fresh, pleated and new. They are purchased specifically for your area so they are always the right length and size, avoiding issues with the wrong reusable curtains being returned from the laundry.

FROM ONLY
\$20.95 + GST
EACH

Reduce annual costs by up to 97%

Contact us for a copy of the independent case study



Inhibits & reduces bacterial growth

Contact us for a copy of the independent laboratory testing results.



Identification Label

Clear date identification allows for improved changeover compliance and infection control.



If you are working for a company or organisation within healthcare, aged care or a hospital, free samples are available. To request a free sample contact Haines® Medical



1.1 PRICES CURRENT FROM 1ST JULY 2019

HAINES® MEDICAL AUSTRALIA. ALWAYS THINKING. ALWAYS THERE.


Specialists in Disposable and Single Patient Use products

26 Heath Street, Lonsdale
South Australia 5160

T 08 8294 5999
F 08 8294 4337

E sales@hainesmedical.com.au
W hainesmedical.com.au

Disposable Curtains



HAINES® DISPOSABLE CURTAINS

Significant cost savings and improved infection control



HAINES® DISPOSABLE CURTAINS - RECYCLABLE

Haines® Disposable Curtains are designed to replace your existing reusable curtains in healthcare facilities such as hospitals, day surgeries, medical centres and other healthcare facilities that require the use of curtains. They offer a cost effective alternative to purchasing and washing reusable curtains and with their proven ability to inhibit bacterial growth, they can improve infection control at your facility. The Haines® Disposable Curtains are:

- Manufactured from 100% polypropylene, non-woven fabric – 120gsm
- Fire retardant (AS 2755.2-1985) – Failed to Burn
- Impregnated with anti-bacterial and anti-mildew chemicals including nanometer silver (curtain and mesh)



Previous studies have demonstrated that pathogens such as Vancomycin resistant enterococci (VRE), methicillin resistant staphylococcus aureus (MRSA) and Clostridium difficile are still present on non-disposable curtains even after one week post laundering (Ohl, M. et al 2011). Even with steam cleaning, the curtains still remained contaminated (White, L.F; Dancer, S.J; Robertson, C; 2007).

REDUCE ANNUAL COSTS BY UP TO 97%

Below are figures sourced from a case study ("Review and Report on Haines® Disposable Curtains – Recyclable") whereby a leading Australian Hospital switched to Haines® Disposable curtains in all wards.

Cost per Curtain	Initial Purchase Price	Wash Cost per Day	Total Cost per Year for 71 Curtains	Total Cost after 5 Years
Haines® Disposable 7.5m x 2.0m Curtain	\$34.95 each*	\$0	\$2,656.20* (includes 5 extra curtain changes)	\$13,281*
Reusable 7.5m x 2.0m Curtain	\$1,035.00 each	\$225 (5 curtains washed per day)	\$81,900	\$483,000

The **cost savings by utilising disposable curtains vs non disposable would be \$469,719*** over a 5 year period. In addition, there is a cost recouped for recycling of \$250 over a 5 year period.

*Note: Cost analysis done as at September 2011.

REDUCE OH&S ISSUES

Reduce OH&S Issues - Haines® full size 7.5m Disposable Curtain weighs 1.6kg whereas comparable reusable curtains can weigh up to 10kg.




Easy to install. Curtains are lightweight and pre-packaged in a plastic sleeve.

RANGE OF COLOURS

Haines® Disposable Curtains are available in a range of standard and special order colours:

6 x STANDARD COLOURS

- Same day dispatch



11 x SPECIAL ORDER COLOURS

- Minimum order quantity is 10 Boxes per Colour
- Delivery time 12 weeks





**TESTING ON
BOTH CURTAINS
& CHEMICALS**

TEST RESULTS

Haines® Medical has undertaken anti-bacterial testing with an Australian-based, NATA approved laboratory. Importantly, the tests were undertaken on the actual Haines® Disposable Curtains, not just the chemical themselves.

Bacteria tested	Test method	90 Mins % Reduction	6 Hours % Reduction	24 Hours % Reduction	6 Month % Reduction*	12 Month % Reduction	24 Month % Reduction	Pass
S. aureus	AATC METHOD 100-2004	> 99.98%	> 99.99%	> 99.99%	> 99.99%	> 99.99%	> 99.98%	✓
MRSA		> 99.93%	> 99.99%	> 99.99%	> 99.99%	> 99.99%	> 99.93%	✓
E Coli				> 99.99%	> 99.99%	> 99.99%	> 99.99%	✓
VRE		> 99.94%	> 99.99%	> 99.99%	> 99.99%	> 99.99%		✓
C.diff		> 99.86%						✓

RECOMMENDED CHANGEOVER PROCEDURE FOR HAINES® DISPOSABLE CURTAINS – RECYCLABLE

Haines® Disposable Curtains with their impregnated anti-bacterial and anti-mildew chemicals, including nanometer silver, have a proven ability to inhibit and reduce bacterial growth. This means that the frequency of curtain change may be less than required by untreated reusable curtains (as per the NHMRC 2010 guidelines as indicated in the right hand column below) thus saving time and money.

The table below is an excerpt from an independent review and report on Haines® Disposable Curtains detailing changeover recommendations. Contact Haines® Medical for a full copy of the review.

Area	Recommended changeover frequency of Haines® Disposable Curtains	For comparative purposes, recommended changeover of Reusable Curtains as per [NHMRC 2010]
Very High Risk • Outbreak of infectious disease	Change upon discharge of the patient/client from bed or room and dispose	Change immediately and upon discharge of the patient/client
High Risk • Intensive care units • High dependency unit • Burns unit, renal units • Day surgery • Operating suites • Emergency wards • Wards for patients/clients with highly transmissible diseases requiring transmission based precautions e.g. MROs (VRE, MRSA, S.aureus, E.Coli) • Patients with Clostridium difficile (C. Diff) or Carbapenem - Resistant Enterobacteriaceae(CRE)	Annually - Change annually or if they are visibly soiled or torn	Monthly - Change monthly or change upon discharge if MRO (Multi-Resistant Organisms) present. If visibly soiled change immediately.
Significant Risk • General wards	Annually - Change annually or if they are visibly soiled or torn	6 Monthly - Change twice per year and immediately if visibly soiled
Low Risk • Rehabilitation • Long term care / nursing homes • Office based practice • Medical centre, radiography	Every 2 years - Change every 2 years or if they are visibly soiled or torn	Annually - Change annually and immediately if visibly soiled

Haines® Disposable Curtains are to be changed more frequently if required based on the healthcare facilities risk assessment or the healthcare facilities policy.

SPECIFICATIONS

- Material: 100% polypropylene non-woven fabric – 120gsm
- Standard eyelets (Option 1) - eyelets are 15cm apart
- Hooks required per Curtain - DC25 x 17; DC45 x 30; DC75 x 50
- Packing: each folded in poly bag
- Tie Back – sewn on so won't fall off
- Identification Label - each curtain is supplied with an easy to read identification label - this is to monitor length of use and provides an easy way to monitor required changes of curtains
- TGA Registered
- Fire Retardant- tested to AS1530.2-1993 and AS2755.2-1985

**6 STANDARD
& 12 SPECIAL
ORDER
COLOURS**

Option 1: Disposable Curtains to fit existing hooks



Summer Blue

Cream

Latte

Grey

Lilac

Paediatric

2.5m x 2m drop

Summer Blue Code: **DC25B**
Cream Code: **DC25C**
Latte Code: **DC25L**
Grey Code: **DC25GY**
Lilac Code: **DC25LILAC**

BOX 15 **\$314.25** (\$20.95 ea)
 \$345.68 INCL GST

Paediatric Code: **DC25PAED**

BOX 15 **\$375** (\$25 ea)
 \$412.50 INCL GST

4.5m x 2m drop

Summer Blue Code: **DC45B**
Cream Code: **DC45C**
Latte Code: **DC45L**
Grey Code: **DC45GY**
Lilac Code: **DC45LILAC**

BOX 8 **\$239.60** (\$29.95 ea)
 \$263.56 INCL GST

Paediatric Code: **DC45PAED**

BOX 8 **\$320** (\$40 ea)
 \$352 INCL GST

7.5m x 2m drop

Summer Blue Code: **DC75B**
Cream Code: **DC75C**
Latte Code: **DC75L**
Grey Code: **DC75GY**
Lilac Code: **DC75LILAC**

BOX 5 **\$199.75** (\$39.95 ea)
 \$219.73 INCL GST

Paediatric Code: **DC75PAED**

BOX 5 **\$325** (\$65 ea)
 \$357.50 INCL GST

Option 2: Hooks

a. Roller Hooks

Replacement curtain track hooks with simple to attach hook for quick curtain fitting and roller wheel for smooth tracking



Code: **DCROLL-50**

PACK 50 **\$35** (\$38.50 INCL GST)

b. "S" Type Hooks

Used to connect Disposable Curtains to existing curtain hooks. Ideal for window applications.



Code: **DCSHOOK-50**


PACK 50 **\$30** (\$33 INCL GST)

Packs of 50 required per BOX of Curtains:

DC25 - 6 packs
 DC45 - 5 packs
 DC75 - 5 packs

Option 3: Disposable Curtains with Mesh

Curtains with 50cm mesh at the top of the curtain for ceiling mounted curtains if compliance with fire regulation is an issue.



Summer Blue

Cream

Latte


4.5m x 2.3m drop

Summer Blue Code: **DC45BM**
Cream Code: **DC45CM**
Latte Code: **DC45LM**
Grey Code: **DC45GYM**
Lilac Code: **DC45LILACM**

BOX 8 **\$300** (\$37.50ea)
 \$330 INCL GST

Paediatric Code: **DC45PAED-M**

BOX 8 **\$370** (\$46.25 ea)
 \$352 INCL GST



Grey

Lilac

Paediatric

7.5m x 2.3m drop

Summer Blue Code: **DC75BM**
Cream Code: **DC75CM**
Latte Code: **DC75LM**
Grey Code: **DC75GYM**
Lilac Code: **DC75LILACM**

BOX 5 **\$325** (\$65ea)
 \$357.50 INCL GST

Paediatric Code: **DC75PAED-M**



BOX 5 **\$360** (\$70 ea)
 \$396 INCL GST

HAINES® MEDICAL AUSTRALIA. ALWAYS THINKING. ALWAYS THERE.

Specialists in Disposable and Single Patient Use products

26 Heath Street, Lonsdale
 South Australia 5160

T 08 8294 5999
F 08 8294 4337

E sales@hainesmedical.com.au
W hainesmedical.com.au  

Disposable Curtains

