

REUSABLE TOURNIQUET

For improved infection control

Haines®
MEDICAL AUSTRALIA

Reusable Tourniquet

FEATURES

- Soft and comfortable band to prevent skin tears and bruising
- Fast and easy application
- Simple click buckle mechanism for securing in place
- Thumb grip for easy strap adjustments
- Auto tension locks helps prevent damage from overtightening
- Light colour enables visualisation of blood or soiling
- Latex free
- Reduces infection risk between patients
- Boxed to allow easy access to fresh tourniquets
- Ensures a firm grip unlike many disposable tourniquets
- 50cm long x 2.5cm wide

Code: **RT-25**

BOX 25 **\$98.75** (\$3.95 ea) GST FREE




Code: **RT-500** (Carton 500 - 20 Boxes of 25)

BOX 500 **\$1500** (\$3 ea) GST FREE



USING A HAINES® REUSABLE TOURNIQUET

- 1 Wrap the Reusable Tourniquet around the patient's arm and click buckle to secure in place
- 2 Pull tight for a firm grip. The auto tension lock helps prevent damage from overtightening
- 3 Press the button to release the buckle for quick and easy removal of the Reusable Tourniquet



Haines® Medical also offers a simple, cost effective **Single Patient Use Tourniquet**



Haines® Medical also offers a simple, cost effective **Single Patient Use Tourniquet:**

Code: **SPUT-25** (Box 25) \$1.99 ea GST FREE

Code: **SPUT-500** (Carton 500 - 20 Boxes of 25) 99c ea GST FREE



If you are working for a company or organisation within healthcare, aged care or a hospital, free samples are available. To request a free sample contact Haines® Medical





PRICES CURRENT FROM 1ST JULY 2018

HAINES® MEDICAL AUSTRALIA. ALWAYS THINKING. ALWAYS THERE.

Specialists in Disposable and Single Patient Use products

26 Heath Street, Lonsdale
South Australia 5160

T 08 8294 5999
F 08 8294 4337

E sales@hainesmedical.com.au
W hainesmedical.com.au  

Tourniquets

