# SALLYSLING<sup>®</sup> GENERAL PURPOSE WITH HEAD SUPPORT

Available with Clip or Loop attachments





## Single patient use

Reduces the risk of cross contamination between patients







Patient comfort Soft polyester material for patient comfort



## Multi-purpose Soft polyester material for patient comfort



# Available in clips or loops

Slings are available with loop or clip attachments to fit hook or loop style lifters

SallySling<sup>®</sup> General Purpose with Head Support slings support the whole body including the head. They are designed for use on one patient during their stay and should not be redistributed to other patients. With built-in head support and padded legs, the SallySling<sup>®</sup> General Purpose with Head Support sling ensures superior patient comfort during transfers. The SallySling<sup>®</sup> General Purpose with Head Support sling can be used for bed, chair and floor lifts.

- Slings are available with loop or clip attachments to fit hook or loop style lifters
- Tested to Australian Standards AS/NZS ISO 10535:2011
- Made of 70% Polyester / 30% Viscose with Polypropylene Straps
- Always ensure you have read the Sling Application Guide before use

# BENEFITS FOR YOUR STAFF AND PATIENTS

- Can be used in pools, bathing and showers
- Soft Poly / Viscose to maintain skin integrity
- Breathable this reduces sweating to help reduce moisture related skin integrity (see air permeable report)
- Absorbent fluids don't run off
- Minimal stretch under load weight
- Can be used during XRAY / MRI



If you are working for a company or organisation within healthcare, aged care or a hospital, free samples are available. To request a free sample contact Haines<sup>®</sup> Medical

HAINES<sup>®</sup> MEDICAL AUSTRALIA. ALWAYS THINKING. ALWAYS THERE. Specialists in Disposable and Single Patient Use products

26 Heath Street, Lonsdale	
South Australia 5160	

T 08 8294 5999F 08 8294 4337

E sales@hainesmedical.com.au W hainesmedical.com.au in f





SallySlings & Transfer Boards

ISO 900' CERTIFIED

# Sallysling® General Purpose Vith Head Support



Available with Clip or Loop attachments

#### **Identification Label**

The sewn-on label has space to write the patient's name to clearly identify the sling as theirs, to reduce any risk of cross-infection.

#### Wash Indicator

The SallySling<sup>®</sup> General Purpose is NOT designed to be washed. If the sling is accidentally washed or



heavily soiled, the wash indicator will change labels so that the care provider is aware the sling not be used further.

#### **Expected Lifetime**

This sling can be used as long as it does not display any signs of wear or damage prior to each use, however we recommend that you replace the sling after 6 months of use

#### **Easy Identification**

Size is indicated by the colour of the binding

#### **Extra Strength**

Unlike some disposable slings, the SallySling® is Polyester Viscose for extra strength

#### Padded legs for comfort

#### Attachment types

South Australia 5160



Loops

	Code	Size	Body	Leg	Width	Binding Colour	Safe Working Load	QTY	Price (GST FREE)
	SLSPU-S	SMALL	77cm	37cm	67cm	RED	200 kg / 440 lbs	BOX <b>10</b>	\$350 (\$35 ea)
Š	SLSPU-M	MEDIUM	85cm	46cm	74cm	YELLOW	200 kg / 440 lbs	BOX 10	\$350 (\$35 ea)
LOOPS	SLSPU-L	LARGE	96cm	56cm	86cm	GREEN	230 kg / 507 lbs	BOX 10	\$350 (\$35 ea)
	SLSPU-XL	X-LARGE	116cm	61cm	109cm	BLUE	230 kg / 507 lbs	BOX <b>10</b>	<b>\$390</b> (\$39 ea)
	SLSPU-XXL	XX-LARGE	127cm	61cm	109cm	ORANGE	400 kg / 881 lbs	BOX <b>5</b>	<b>\$250</b> (\$50 ea)
	SLSPUCLIP-S	SMALL	81cm	33cm	78cm	RED	200 kg / 440 lbs	BOX <b>10</b>	\$350 (\$35 ea)
10	SLSPUCLIP-M	MEDIUM	90cm	38cm	90cm	YELLOW	200 kg / 440 lbs	BOX <b>10</b>	\$350 (\$35 ea)
CLIPS	SLSPUCLIP-L	LARGE	102cm	50cm	102cm	GREEN	230 kg / 507 lbs	BOX <b>10</b>	<b>\$350</b> (\$35 ea)
	SLSPUCLIP-XL	X-LARGE	117cm	50cm	108cm	BLUE	230 kg / 507 lbs	BOX <b>10</b>	<b>\$390</b> (\$39 ea)
	SLSPUCLIP-XXL	XX-LARGE	132cm	60cm	124cm	ORANGE	400 kg / 881 lbs	BOX 5	<b>\$250</b> (\$50 ea)


Always ensure you have read the sling application g

Note: The XXL Sling is a 2 layer white polyester mesh construction. Shape and

#### HAINES® MEDICAL AUSTRALIA. ALWAYS THINKING. ALWAYS THERE. Specialists in Disposable and Single Patient Use products 26 Heath Street, Lonsdale **T** 08 8294 5999 E sales@hainesmedical.com.au

**F** 08 8294 4337

		4 = - ((				
200 kg / 440 lbs	BOX <b>10</b>	<b>\$350</b> (\$35 ea)				
200 kg / 440 lbs	BOX 10	\$350 (\$35 ea)				
230 kg / 507 lbs	BOX <b>10</b>	\$350 (\$35 ea)				
230 kg / 507 lbs	BOX <b>10</b>	<b>\$390</b> (\$39 ea)				
400 kg / 881 lbs	BOX <b>5</b>	<b>\$250</b> (\$50 ea)				
guide before use						
d features vary to the above information						
-						



W hainesmedical.com.au in f SallySlings & Transfer Boards