
HANDICARE'S PRISM Sling Selection






COMFORT RECLINE SLING


The Comfort Recline Sling is designed for the amputee client but it is also suitable for clients who have to remain in the sling for extended periods. The spacer material is fully breathable and stretches to mold to the shape of a chair. The sling features elasti-cated sections which allow the straps to be folded safely out of the way when the client is seated.



The Universal Sling is the most commonly supplied sling for general transfer purposes. The universal is easy to fit yet is very sup-portive. It can be used in a divided leg or crossover position.








HAMMOCK SLING

The Hammock Sling provides maximum support to the sacral area and is used in cases where the universal sling is unsuitable because the individual has poor body tone. The sling can be used in a divided leg or crossover position; this also provides good support for amputees.

DUAL ACCESS SLING

The Dual Access Sling is designed to be used with either a ceiling lift or a stand aid. It gives access to be able to remove and refit lower garments for toileting transfers while the client is in the sling.





CUSTOM SLING

Although we offer over 100 standard slings, everyone is unique and their needs are different. That is why we will customize any of our standard slings to help suit your needs better.

You can modify an existing sling design, like adding an extra strap, or we can start right from scratch.

Before creating a custom sling, our engineering department reviews the design to ensure that there will be no loss of integrity with the changes.

Once the custom sling has been made by our in-house sling manufacturing, we weight test each custom sling to ensure its strength.

Common customizations include adding: loops, straps and padding.

CLIENT SPECIFIC SLING

With Infection Control becoming an area of greater focus and attention, Client Specific Slings are designed to be used by one individual and disposed of when the patient/resident no longer requires or the sling is soiled.

The Client Specific slings are available in a variety of styles and sizes. These slings are made of a strong, non-woven material that conforms to the body to improve comfort ad safety. And because these slings are intended for only one client, they help prevent infections from spreading as well as eliminating laundry costs.



Sling Sizing Chart

Size	Height	Weight
Junior	<4'	<110 lbs
Small	4'-5'6"	95 lbs-150 lbs
Medium	5'-6'	125 lbs-250 lbs
Large	6'-7'	250 lbs-400 lbs
X-Large	6'-7'	400 lbs-480 lbs
XX-Large	> 6.5'	>480 lbs



WE'RE THE COMPANY THAT BRINGS YOU A Wide Selection of Slings

EXCEPTIONAL FIT AND DURABILITY



SLINGS TO ACCOMMODATE VIRTUALLY ANY LIFT, POSITIONING AND TRANSFER REQUIREMENT

Handicare's Prism sling line is recognized by healthcare professionals throughout North America for exceptional fit and durability. Our range of sling styles and sizes will accommodate virtually any lift, positioning and transfer requirement. And, if our standard selection does not provide a sling solution, we will custom design a sling system to your specifications.

Prism slings are manufactured from durable nylon mesh, polyester fabrics and quilted padding to provide contoured support, comfort and safety. As well, these long-wearing, quick-drying slings contain no metal parts, making them convenient to launder and easy to return to service, without the need to locate mislaid inserts.




10888 Metro Court, St. Louis, MO 63043 | 866.891.6502 | www.handicareusa.com 485 Millway Avenue, Unit 2, Concord, ON L4K 3V4 | 877.304.5438 | www.handicare.ca

HAMMOCK SLING

Acknowledged by professionals as the most comfortable and supportive sling for both home and institutional use, the versatile Hammock Sling combines toileting access with full body and head support. The sling comes with leg and head straps and features loops to allow for various sitting positions. The Hammock is available in either regular mesh or with quilted material.









HYGIENE SLING

The Hygiene Sling was specifically designed for toileting and hygiene functions. Access is unrestricted for most users, making this sling ideal for independent transfer and toileting. The Hygiene Sling features quilted material and comes with leg and shoulder straps with loops for various seating positions. A wide velcro strap across the chest provides additional upper body support.





DELUXE WALKING SLING


The Deluxe Walking Sling is the new generation of slings for support during ambulation. This sling features all of the same features of the Walking Sling, but offers additional adjustment around the torso.



UNIVERSAL SLING

The Universal Sling is easy-to-use and easy-to-fit. It is a multi-purpose sling that provides toileting access as well as good back and thigh support. The sling comes with leg, hip and shoulder straps and features loops to accommodate various seating positions. Available in either regular mesh or with quilted material, the Universal Sling can also be supplied with added head and neck support.






STAND-AID SLING

Designed by our expert research team, the Stand-Aid sling will comfortably wrap around the patients waist and with limited pressure will help raise the person to a standing position. The Stand-Aid Sling is padded for extra comfort and comes with Prism's signature 3 loop design, which allows for greater options for comfort when performing the lift.







REHAB WALKING SLING

Prism's Rehab Walking Sling provides a natural sit to stand movement with non-slip fabric to keep the patient stationary while moving. The Rehab Walking Sling is functional from both bedside and chair surface and encourages pivot transfers, ambulation and toileting.



WALKING SLING

The Walking Sling was developed with the help of physiotherapists. This supportive sling provides lateral stability during gait training. The Walking Sling offers upper body support and freedom of movement to individuals performing exercises in a rehabilitation setting.



STANDING SLING


The Standing Sling is a dual purpose sling affording a natural sit to stand movement. Functional from both bedside and chair surface, this sling encourages pivot transfers and ambulation.



POSITIONING SLING

The Positioning Sling is intended to be permanently in place as part of the bedding. The sling has six pairs of straps and features loops to support the head in various positions. Repositioning and turning, done manually under time pressure, poses a significant risk of injury to caregivers. The Positioning Sling is designed to address that risk. The sling is also used for transfers to other horizontal surfaces.









TRI-TURNER SLING

Ideal for turning a client into a lateral, resting position, the Tri-Turner sling also has a sewn-in velcro loop & hook to attach the sling to the bed, providing a more comfortable turn for the client.



Ideal for supporting limbs to facilitate change of dressing on arms, hands, legs, feet or to elevate and hold limbs during surgery.

