



STEARIC ACID GUN KIT INSTRUCTIONS

Roger George's Stearic Acid Gun Kit comes with:

- 1 - mold with adjustable depth
- 1 - small jar of stearic acid
- 1 - stainless steel spoon
- 1 - 2-oz eyedropper bottle of Roger George blood
- 3 - razor blades
- 1 - blow gun

To make a hit:

- Set desired depth on your mold by moving the shaft up or down. When correct depth is achieved, tighten knob to lock shaft in place
 - Max depth: 1/2 inch
- Heat a spoonful of stearic acid until it melts
 - Heat from the bottom center of the spoon to keep the acid from splattering
- Once the acid is melted, pour it into the mold
- Count to 12, then pour the residual acid back into the spoon
- Use razor blade to clean excess acid off surface of mold
- Fill the remaining void with drops of blood
 - Do not overfill - leave room for stearic acid to cap
 - 9-10 drops is the max you should use
- Pour the remaining liquified stearic acid into mold to create a cap over the blood
- Allow a minute or so for acid to set, then scrape excess off mold with razor
- Loosen knob enough for the hit to move and gently ease the shaft up until the hit is free of the mold
- Push the hit off the top of the shaft being careful not to squeeze it
- Load the hit into the blow gun from the back near the mouth piece
- Fire away
 - You'll probably need to try a few test shots to get the hang of it
 - Get close to the person you're shooting at for best effect

The stearic acid gun is a non-pyrotechnic product and does not require a license.

For a demonstration please check out the video on our website:

<https://rogergeorge.com/StearicAcidGunVideo>