

## STRENGTHEN & REPAIR NAILS

# IBX Treatment

IBX is a strengthening system that toughens, repairs, and protects natural nails. This effective treatment is used to improve nail plate integrity under any gel polish or enhancement service. It can also be used to target natural nail imperfections and weaknesses that prevent them from growing long.



Treatment time:

**20 min.**



Suggested price:

**\$25**



Cost per treatment

**35¢**



### SAMPLE MENU DESCRIPTION:

The IBX Treatment starts with IBX Repair to address peeling, splits, and free edge delamination for damage control and protection under gel polish. Repair works like "double-sided sticky tape" to seal damaged areas together prior to the next step, the application of IBX Strengthen.

IBX Strengthen will create nail strength from within, under gel polish. It's a penetrating toughening agent that works inside the nail plate.



# What you'll need:



**IBX STRENGTHEN, 0.5 OZ**  
BX4010

Repair dull, cracked nails with IBX Strengthen. IBX is the first of its kind, working inside the nail, not on top. It penetrates into the upper layers of the nail and cross-links inside the nail itself creating an interpenetrating polymer network.



**IBX REPAIR, 0.33 OZ**

BX4015



**CINNATIZE, 32 OZ**

N2670



**INTRINSICS LINT FREE NAIL WIPES, 200 CT**

CX122T



**CUCCIO BLACK MINI LED CURING LAMP**

N3534



**IKONNA 5" WHITE NAIL FILES, 100/180 GRIT, 50 PC**

NB225



**CHROME PLATED TOE NAIL CLIPPER, STRAIGHT**

NB135



# Protocol

## IBX Repair



1

### APPLY

After you have removed polish, filed, and shaped nails, you will cleanse the nail with a lint-free pad and Cinnatize. Gently shake IBX Repair treatment, you will need to apply IBX Repair twice.



2

### HEAT

Gentle heat for 1 minute. You may use a nail polish dryer set to warm or a table lamp with a 75-100 watt bulb.



3

### BLOT

Lightly blot the excess from the surface if needed.



4

### CURE

Cure with an LED or UV Lamp for the same amount of time your lamp cures your Gel Polish color coat.



5

### CLEANSE

Cleanse with your preferred cleanser. We recommend Cinnatize.



# Protocol

## IBX Strengthen



1

### APPLY

Gently shake IBX Strengthen treatment and apply a thin layer over the nail plate.



2

### HEAT

Gentle heat for 4 minutes. You may use a nail polish dryer set to warm or a table lamp with a 75-100 watt bulb.



3

### BLOT

Lightly blot the excess from the surface if needed.



4

### CURE

Cure with an LED or UV Lamp for the same amount of time your lamp cures your Gel Polish color coat.



5

### CLEANSE

Cleanse with your preferred cleanser. We recommend Cinnatize.

\*Repeat the 5 steps above for a second application of IBX.