

Prescription Approval Form

Please complete this page and return it via fax (647-905-7119) or email (info@vetmedsolutions.com). If this form is received Monday-Friday between 9am-5pm EST, we will contact the pet owner within 2 hours.

PATIENT INFORMATION

Owner Name: _____

Patient Name: _____ Age: _____

Owner Phone: _____

Species/Breed: _____

Owner Email: _____

Weight: _____ Sex: _____

Owner Address: _____

Relevant Conditions: _____

PRESCRIPTION INFORMATION

Portable Oxygen Kit

Portable Oxygen (99%), Flow by Mask

Choose a Kit:

- | | |
|---------------------------------------------------|--------------------------------------------------|
| <input type="checkbox"/> Dog, 0-11 lb (0-5 kg) | <input type="checkbox"/> Cat, 0-11 lb (0-5 kg) |
| <input type="checkbox"/> Dog, 11-22 lb (6-10 kg) | <input type="checkbox"/> Cat, 11-22 lb (6-10 kg) |
| <input type="checkbox"/> Dog, 22-33 lb (11-15 kg) | <input type="checkbox"/> Custom Kit: _____ |
| <input type="checkbox"/> Dog, 33+lb (16+ kg) | _____ |

Instructions for use: _____

Extended Oxygen Therapy

Oxygen (>95%), Flow by Mask or Cage

PLEASE INDICATE APPROVAL:

- APPROVED NOT APPROVED

Includes an Oxygen Concentrator (>95% purity @ 1- 10 LPM), an Oxygen Cage, and a Pet Oxygen Mask. Flow rate depends on the delivery method. Pet owners will be informed using guidelines on Page 2. Veterinarians may write further instructions below:

HOSPITAL INFORMATION

Hospital Name: _____

Contact Number: _____

License: _____ Province: _____

Contact Email: _____

Authorizing Doctor: _____

(print)

(sign)

Fax to 647-905-7119 or email to info@vetmedsolutions.com. Questions? Call (888) 976-8696 today!

For more information about prescribing Pawprint Oxygen, please visit www.vetmedsolutions.ca

Prescription Guidelines

Veterinarians: You are responsible for determining the oxygen flow rate based on the delivery method, animal weight, and desired inspired concentration. While Pawprint Oxygen provides guidelines (see chart below), it is the attending veterinarian's responsibility to determine dosage and oxygen flow rate.

Portable Oxygen Kit

A portable supply of oxygen for use at home and during transport. Each Kit includes Pawprint Oxygen Canisters, Regulator, PureVent Mask and Tubing.

Standard Kits:

Kit	Mask	Flow Rate	Duration
Cat 0-11lb	Small	0.5 LPM	90 min
Cat 11+ lb	Small	0.1 LPM	90 min
Dog 0-11 lb	Medium	0.5 LPM	90 min
Dog 11-22 lb	Medium	0.1 LPM	90 min
Dog 22-33 lb	Large	0.2 LPM	45 min
Dog 33+ lb	Large	0.3 LPM	30 min

The **Flow Rate** for standard kits is based on a dosage of 100 ml/kg/min. The Regulator is factory-set to the prescribed flow rate and cannot be changed by the client.

The **Duration** of the kit can be increased by purchasing additional Pawprint Oxygen Canisters, which contain 15 liters of oxygen (99%) each.

The **Mask** refers to the PureVent Pet Oxygen Mask, which is included in all kits.

Custom Kits: Write in the flow rate (0.5,1.0, 2.0, or 3.0 LPM). The maximum flow rate is 3.0 LPM. If you wish to prescribe a higher flow rate, please consider an Oxygen Concentrator. For custom kits, the correct Mask size will be provided by our experts based on the pet.

Extended Oxygen Therapy

Oxygen Concentrators may be used for pets that require more than one hour of therapy per day or pets requiring higher flow rates than can be achieved with Oxygen Canisters.

Concentrators may be used in tandem with a Mask or a Cage.

The Concentrator is available in two models:

- 5 Liter: Adjustable flow rate from 1-5 LPM
- 10 Liter: Adjustable flow rate from 2-10 LPM

Oxygen Tents are available in the following sizes:

- Small (10" x 8" x 8"): Rabbits / exotics
- Medium (20" x 13" x 14"): Cats / small dogs.
- Large (31" x 21" x 24"): Middle sized dogs.

To achieve sufficient oxygen saturation levels::

- Use a 5 Liter Concentrator to supply oxygen to a Small or Medium Tent.
- Use a 10 Liter Concentrator to supply oxygen to a Large Tent

Note: Pawprint Oxygen Canisters cannot be used to supply Oxygen Tent due to flow rate requirements.

Ordering Process

VetMed Solutions is open Monday-Friday 9am-5pm EST. If you are submitting this form within normal business hours, pet owners can expect a call from our team within 2 hours to complete their order. If you are submitting this form outside of normal business hours, pet owners can expect a call the next business day. If this is urgent or requires special instructions, please call us at 888-976-8696. If you would like to dispense Pawprint Oxygen directly from your practice, please call or visit vetmedsolutions.ca