

For accuracy, please ensure that PDF is printed at 100%



We make our rings in 2 sizes :  
MED and LGE

- MED - 19 mm Inner diameter.  
Inner circumference - 59.5  
UK size S. US size 9.
- LGE - 20.6mm inner diameter.  
Inner circumference - 64.6  
UK size W. US size 11


To find your size, choose from the following 3 methods.

### THE STRING TEST

- Cut a small piece of string, thread or ribbon from a spool.
- **MEASURE:** Wrap it snugly around your finger for an accurate measurement.
- **DETERMINE:** Select the most suitable size according to inner circumference measurements on the left.
- **Tip:** Be careful not to wrap the string too tight. Aim for a comfortable, snug fit.
- **Tip:** For the most accuracy, have someone help you measure.
- **CHECK:** Double check for accuracy.

### THE PAPER TEST

- Print and carefully cut out our paper ring sizer below.
- Slip the pointed end through the slit and place around your finger.
- **MEASURE:** Pull the pointed end so it slides to a snug fit.
- **Tip:** Be careful not to wrap the paper too tight. Aim for a comfortable, snug fit.
- **Tip:** For the most accuracy have someone help you measure.
- **CHECK:** Double check for accuracy.



### THE PLACEMENT TEST

- Print ring size circles.
- **MEASURE:** Place a ring that you already wear on top of each circle until you find a perfect match.
- **DETERMINE :** The circle line should match the inner circumference of your ring exactly
- **CHECK:** Double check for accuracy.

